

---

**0090005      Acetone, Quantitative****AC****Performed:** Sun-Sat  
**Reported:** 1-3 days

**Specimen Required:** Patient Prep: For medical purposes only. Timing of specimen collection: Dependent on time of exposure, test upon presentation to hospital.  
Collect: Plain Red. Also acceptable: Lavender (EDTA), Pink (K<sub>2</sub>EDTA), Green (Sodium or Lithium Heparin), or Gray (Potassium Oxalate/Sodium Fluoride).  
Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 0.5 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.2 mL) **Cap tube tightly to minimize alcohol loss.**  
Storage/Transport Temperature: Refrigerated.  
Unacceptable Conditions: Whole blood. Plasma Separator Tubes (PST), Serum Separator Tubes (SST).  
Stability (collection to initiation of testing): After separation from cells: Ambient: 5 days; Refrigerated: 1 week; Frozen: 1 month