
2007949 Paliperidone, Serum or Plasma

PALIPERID

Performed: Wed, Sat
Reported: 1-5 days

Specimen Required: Patient Prep: Pre-dose (trough) draw - At steady state concentration.

Collect: Lavender (EDTA), Pink (K₂EDTA), or Plain Red.

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.5 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Gel separator tubes, light blue (citrate), or yellow (SPS or ACD solution).

Stability (collection to initiation of testing): Ambient: 2 weeks; Refrigerated: 2 weeks; Frozen: 2 months

HOTLINE NOTE: Remove information found in the Specimen Required Remarks field. There is also a component change associated with this test.

Remove component 2011545, Paliperidone Dose

Remove component 2011546, Paliperidone Dose Frequency

Remove component 2011547, Paliperidone Route

Remove component 2011548, Paliperidone Type of Draw