

Quarterly HOTLINE: Effective February 20, 2018

New Test **3000142** **Alpha Fetoprotein (Amniotic Fluid) with Reflex to Acetylcholinesterase and Fetal Hemoglobin** **AF AFP**



Patient History for Prenatal Cytogenetics

Methodology: Quantitative Chemiluminescent Immunoassay/Electrophoresis
Performed: Sun-Sat
Reported: 3-4 days
 Reflex: 3-11 days

Specimen Required: Patient Prep: Amniocentesis. Specimen must be drawn between 13 weeks, 0 days and 36 weeks, 6 days gestation.
Collect: Amniotic fluid.
Specimen Preparation: Transport 2.5 mL amniotic fluid. (Min: 1.5 mL)
Storage/Transport Temperature: Room temperature.
Remarks: **Submit with Order:** Gestational age at time of collection or estimated due date.
Unacceptable Conditions: Specimens contaminated with fetal blood.
Stability (collection to initiation of testing): Ambient: 1 month; Refrigerated: 3 months; Frozen: 3 months

Reference Interval:

Test Number	Components	Reference Interval
	AFP, Amniotic Fluid	By report Ranges are based upon the weeks of gestation.
2006848	Acetylcholinesterase and Fetal Hemoglobin, Amniotic Fluid	Acetylcholinesterase: Negative Fetal Hemoglobin: Negative
	Multiple of Median	1.99 or less

Interpretive Data: Refer to report.

Note: Information must include weeks of gestation. If the AFP (amniotic fluid) is elevated, then Acetylcholinesterase will be added. Additional charges apply. Acetylcholinesterase testing requires an additional 3-11 days to be reported.

CPT Code(s): 82106; if reflexed, add 82013 and 83033

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.