

New Test
Available Now

3000486

Aldosterone Right Adrenal Vein

ALDO RAV



Additional Technical Information

Methodology: Quantitative Chemiluminescent Immunoassay
Performed: Sun-Sat
Reported: Within 24 hours

Specimen Required: Patient Prep: Adrenal venous sampling procedure is required.

Collect: Serum Separator Tube (SST) or Plain Red.

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)

Storage/Transport Temperature: Frozen.

Unacceptable Conditions: EDTA plasma.

Stability (collection to initiation of testing): After separation from cells: Ambient: 8 hours; Refrigerated: 5 days; Frozen: 1 month

Note: Refer to the Additional Technical Information for Endocrine Society recommendations for patient preparation, specimen collection, medications for hypertension control during confirmatory testing for primary aldosteronism, and factors that may lead to false-positive or false-negative aldosterone-renin ratio (ARR) results.

CPT Code(s): 82088

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.