

# BRAF V600E Mutation in Hairy Cell Leukemia

## Indications for Ordering

---

- Confirm diagnosis of hairy cell leukemia (HCL)
- Monitor tumor burden

## Test Description

---

- Genomic DNA is extracted
- Polymerase chain reaction (PCR) amplification of fragment spanning the *BRAF* V600 codon with allele-specific primers for the wild type and the *BRAF* V600E mutant allele
- Quantitation using hydrolysis probe
- Relative percentages of the wild type of *BRAF* V600 and V600E mutant alleles are calculated using a heterozygous calibrator plasmid

## Tests to Consider

---

### Primary test

[BRAF V600E Mutation Detection in Hairy Cell Leukemia by Real-Time PCR, Quantitative 2007132](#)

- Diagnosis/monitoring of HCL

### Related test

[Leukemia/Lymphoma Phenotyping by Flow Cytometry 2008003](#)

- Initial testing to establish tumor lineage

## Disease Overview

---

**Prevalence** – rare lymphoproliferative disorder

### Diagnostic issues

*BRAF* V600E is a reliable molecular marker to confirm diagnosis of HCL

- Mutations detected in nearly all cases of HCL but rarely in other lymphoproliferative disorders (Tiacchi E, 2011)

### Treatment issues

Quantitation of allele burden allows monitoring of response to therapy

## Genetics

---

**Gene** – *BRAF*

### Structure/function

- *BRAF* protein kinase acts in the RAS/mitogen-activated protein kinase-signaling pathway
- Major role in cell proliferation, survival, and neoplastic transformation

### Mutations

Most mutations occur at codon V600

- Mutation results in V600E change

## Test Interpretation

---

**Analytical sensitivity** – 0.2% mutant allele

### Results

- Positive – *BRAF* V600E allele detected and quantified
- Weakly positive, nonquantifiable – *BRAF* V600E mutation detected at 0.2-0.4% mutant allele

### Limitations

Limit of detection is 0.2% mutant allele