

APC- and MUTYH-Associated Polyposis Panel, Sequencing and Deletion/Duplication

Last Literature Review: November 2021 Last Update: June 2023

Familial adenomatous polyposis (FAP) is an *APC*-associated polyposis condition caused by pathogenic variants in the *APC* gene resulting in the development of hundreds to thousands of adenomatous colonic polyps beginning in early adolescence. The lifetime risk for colorectal cancer (CRC) in individuals with FAP is 100%. Additional symptoms may include dental anomalies, polyps of the gastric fundus and duodenum, and congenital hypertrophy of the retinal pigment epithelium (CHRPE). Pathogenic *APC* variants may also cause other related polyposis conditions, including attenuated FAP (AFAP) or gastric adenocarcinoma and proximal polyposis of the stomach (GAPPS).

MUTYH-associated polyposis (MAP), caused by biallelic pathogenic variants in the *MUTYH* gene, can result in the development of colon polyps that are less numerous (typically 10-100s) and is often diagnosed later in life. Genetic testing may be used to assess individuals at risk for FAP, other *APC*-associated polyposis conditions, or MAP due to a suggestive personal or family history.

Disease Overview

Associated Disorders

Disorder	Polyp and Cancer Characteristics	Age of Onset	Other Symptoms
FAP	<p>Hundreds to thousands of adenomatous colonic polyps may develop</p> <p>Polyps may develop in gastric fundus and duodenum</p>	<p>Symptom onset at 7-36 years (generally early adolescence)</p> <p>Mean age of CRC diagnosis is 39 years in untreated individuals</p> <p>Without preventive colectomy, all individuals will develop CRC</p>	<p>Dental anomalies</p> <p>CHRPE</p> <p>Osteomas, soft tissue tumors, desmoid tumors</p>
AFAP	<p>Fewer colonic polyps than FAP (10-100s, with an average of 30)</p> <p>More proximally located polyps</p>	<p>Cancer occurs later than in FAP</p>	<p>Extracolonic manifestations are variable</p>
GAPPS (associated with pathogenic variants in promoter 1B of the <i>APC</i> gene)	<p>Gastric fundic gland polyposis and increased risk for gastric cancer</p> <p>Limited colonic involvement</p>	<p>Unknown</p>	<p>Unknown</p>
MAP (biallelic pathogenic <i>MUTYH</i> variants)	<p>10 to a few hundred colonic adenomatous polyps</p>	<p>Third decade or later</p>	<p>Unknown</p>

Source: Jasperson, 2017¹; Nielsen, 2021²

Featured ARUP Testing

To compare directly to other hereditary cancer panels offered by ARUP Laboratories, see the [ARUP Hereditary Cancer Panel Comparison](#) table.

APC- and MUTYH-Associated Polyposis Panel, Sequencing and Deletion/Duplication 3004407

Method: Massively Parallel Sequencing

- Preferred diagnostic or predictive test for APC-associated polyposis conditions (FAP, AFAP, GAPPS, and MAP).
- Testing minors for adult-onset conditions is not recommended; testing will not be performed in minors without prior approval. For additional information, please contact an ARUP genetic counselor at 800-242-2787 ext. 2141.
- If a familial sequence variant has been previously identified, targeted sequencing for that variant may be appropriate; refer to the [Laboratory Test Directory](#) for additional information.

Genetics

Genes

APC (NM_000038, NM_001127511 Exon 1b only) and *MUTYH* (NM_001128425)

For more information, see the [Genes Tested](#) table

Penetrance

Classic FAP: 100% incidence of CRC in untreated individuals¹

MAP:

- 43-63% incidence of CRC by age 60²
- 80-90% lifetime risk of CRC²

Etiology

FAP is estimated to account for 0.5% of CRC cases.¹

MAP is estimated to account for 0.7% of all CRC cases.²

Prevalence

FAP: Approximately 1 in 6,850 to 1 in 31,250 individuals have FAP.¹

MAP: Approximately 1 in 20,000 to 1 in 60,000 individuals have MAP²; approximately 1% of White individuals are predicted to carry a single pathogenic *MUTYH* variant.²

Inheritance

APC is autosomal dominant.¹

MUTYH is autosomal recessive.²

De Novo Variants

APC : 20-25% of individuals with FAP have a de novo pathogenic variant.³

- 20% of individuals with apparent de novo variants in *APC* have somatic mosaicism.¹

Test Description

See the [Genes Tested](#) table for genes included in the panel.

Methodology

This test is performed using the following sequence of steps:

- Selected genomic regions, primarily coding exons and exon-intron boundaries from the targeted genes are isolated from extracted genomic DNA using a probe-based hybrid capture enrichment workflow.
- Enriched DNA is sequenced by massively parallel sequencing (MPS; also known as NGS) followed by paired-end read alignment and variant calling using a custom bioinformatics pipeline.
- Sanger sequencing is performed as necessary to fill in regions of low coverage and, in certain situations, to confirm variant calls.
- The pipeline includes an algorithm for detection of large (single exon-level or larger) deletions and duplications.
- Large deletion/duplication calls made using MPS are confirmed by an orthogonal exon-level microarray when sample quality and technical conditions allow.

Clinical Sensitivity

Classic FAP: approximately 93%⁴

- ≤90% of *APC* pathogenic variants detected by sequencing¹
- Approximately 8-12% of *APC* pathogenic variants detected by deletion/duplication testing^{5,6}

Attenuated FAP: <30%⁷

GAPPS: unknown¹

MAP:

- Approximately 99% of pathogenic *MUTYH* variants detected by full gene sequencing²
- <1% of pathogenic *MUTYH* variants detected by deletion/duplication analysis²

Analytical Sensitivity

Variant Class	Analytical Sensitivity (PPA) Estimate ^a (%) and 95% Credibility Region (%)	Analytical Specificity (NPA) (%)
SNVs	>99 (96.9-99.4)	>99.9
Deletions 1-10 bp ^b	93.8 (84.3-98.2)	>99.9
Insertions 1-10 bp ^b	94.8 (86.8-98.5)	>99.9
Exon-level ^c Deletions	97.8 (90.3-99.8) [2 exons or larger] 62.5 (38.3-82.6) [single exon]	>99.9
Exon-level ^c Duplications	83.3 (56.4-96.4) [3 exons or larger]	>99.9

^aGenes included on this test are a subset of a larger methods-based validation from which the PPA values are derived. These values do not apply to testing performed by multiplex ligation-dependent probe amplification (MLPA).

^bVariants greater than 10 bp may be detected, but the analytical sensitivity may be reduced.

^cIn most cases, a single exon deletion or duplication is less than 450 bp and 3 exons span a genomic region larger than 700 bp.

bp, base pairs; NPA, negative percent agreement; PPA, positive percent agreement; SNVs, single nucleotide variants

Testing Strategy

Contraindications for Ordering

- Should not be ordered to detect somatic variants associated with malignancy as sensitivity for mosaic variants is low with methodology used for germline assays
- Individuals with hematologic malignancy and/or a previous allogeneic bone marrow transplant should not undergo molecular genetic testing on peripheral blood specimen.
 - Testing of cultured fibroblasts is required for accurate interpretation of test results.

Results

Result	Variant(s) Detected	Clinical Significance
Positive	Single pathogenic variant in <i>APC</i> gene	Predictive of FAP or other <i>APC</i> -associated polyposis condition
	Two pathogenic variants in <i>MUTYH</i> gene on opposite chromosomes	Predictive of MAP
	Single pathogenic <i>MUTYH</i> variant	Individual is a carrier of MAP and may be affected if another unidentified pathogenic <i>MUTYH</i> variant is present on the opposite chromosome Possible increased risk for cancer has been associated with a single pathogenic <i>MUTYH</i> variant, but is not well defined

Result	Variant(s) Detected	Clinical Significance
Negative	No pathogenic variants in either <i>APC</i> or <i>MUTYH</i> gene	Does not rule out FAP, other <i>APC</i> -associated polyposis conditions, or MAP
Inconclusive	Variant of uncertain clinical significance detected	Uncertain

Limitations

- A negative result does not exclude a diagnosis of FAP or MAP.
- Diagnostic errors can occur due to rare sequence variations.
- Interpretation of this test result may be impacted if this patient has had an allogenic stem cell transplantation.
- The following will not be evaluated:
 - Variants outside the coding regions and intron-exon boundaries of targeted genes
 - Regulatory region and deep intronic variants
 - Breakpoints of large deletions/duplications
- The following may not be detected:
 - Deletions/duplications/insertions of any size by massively parallel sequencing
 - Large duplications less than 3 exons in size
 - Noncoding transcripts
 - Some variants due to technical limitations in the presence of pseudogenes, repetitive, or homologous regions
 - Low-level somatic variants

Genes Tested

To compare directly to other hereditary cancer panels offered by ARUP Laboratories, see the [ARUP Hereditary Cancer Panel Comparison table](#).

Gene	MIM Number	Disorder/Associated Cancer(s)/Tumor(s)	Inheritance
<i>APC</i>	611731	FAP AFAP GAPPS Associated cancer(s)/tumor(s): colorectal adenomas and cancer, duodenal adenomas and cancer, fundic gland polyps, osteomas, thyroid, pancreas, and others	AD
<i>MUTYH</i>	604933	Associated cancer(s)/tumor(s): breast, ^a colorectal ^a	AD
		MAP Associated cancer(s)/tumor(s): colorectal adenomas and cancer, duodenal adenomas and cancer	AR

^aAssociation is suggested but not well-established at this time.

AD, autosomal dominant; AR, autosomal recessive

References

1. Jasperson KW, Patel SG, Ahnen DJ. [APC-associated polyposis conditions](#). In: Adam MP, Ardinger HH, Pagon RA, et al, eds. *GeneReviews*. University of Washington, Seattle. Last update Feb 2017; accessed Aug 2021.
2. Nielsen M, Infante E, Brand R. [MUTYH polyposis](#). In: Adam MP, Ardinger HH, Pagon RA, et al, eds. *GeneReviews*. University of Washington, Seattle. Last update May 2021; accessed Aug 2021.
3. Bisgaard ML, Fenger K, Bülow S, et al. [Familial adenomatous polyposis \(FAP\): frequency, penetrance, and mutation rate](#). *Hum Mutat*. 1994;3(2):121-125.
4. Lagarde A, Rouleau E, Ferrari A, et al. [Germline APC mutation spectrum derived from 863 genomic variations identified through a 15-year medical genetics service to French patients with FAP](#). *J Med Genet*. 2010;47(10):721-722.
5. Aretz S, Stienen D, Uhlhaas S, et al. [Large submicroscopic genomic APC deletions are a common cause of typical familial adenomatous polyposis](#). *J Med Genet*. 2005;42(2):185-192.
6. Bunyan DJ, Eccles DM, Sillibourne J, et al. [Dosage analysis of cancer predisposition genes by multiplex ligation-dependent probe amplification](#). *Br J Cancer*. 2004;91(6):1155-1159.
7. Lefevre JH, Rodrigue CM, Mourra N, et al. [Implication of MYH in colorectal polyposis](#). *Ann Surg*. 2006;244(6):874-880.

Additional Resources

Li J, Woods SL, Healey S, et al. [Point mutations in exon 1B of APC reveal gastric adenocarcinoma and proximal polyposis of the stomach as a familial adenomatous polyposis variant.](#) *Am J Hum Genet.* 2016;98(5):830-842.

Lubbe SJ, Di Bernardo MChiara, Chandler IP, et al. [Clinical implications of the colorectal cancer risk associated with MUTYH mutation.](#) *J Clin Oncol.* 2009;27(24):3975-3980.

National Comprehensive Cancer Network. [NCCN Clinical Practice Guidelines in Oncology: genetic/familial high-risk assessment—colorectal.](#) Version 1.2021. Updated May 2021; accessed Nov 2021.

Related Information

[Colorectal \(Colon\) Cancer](#)
[Hereditary Gastrointestinal Cancer Panels](#)

ARUP Laboratories is a nonprofit enterprise of the University of Utah and its Department of Pathology. 500 Chipeta Way, Salt Lake City, UT 84108
(800) 522-2787 | (801) 583-2787 | aruplab.com | arupconsult.com