

# Myeloid Malignancies Mutation and Copy Number Variation Panel by Next Generation Sequencing

Last Literature Review: August 2023    Last Update: March 2026

Myeloid malignancies are clonal disorders of hematopoietic stem and progenitor cells that include myelodysplastic syndromes (MDSs), myeloproliferative neoplasms (MPNs), myelodysplastic/myeloproliferative neoplasms (MDS/MPNs), acute myeloid leukemia (AML), and others. Recent studies have identified recurrently mutated genes with diagnostic, prognostic, and therapeutic impact in myeloid malignancies. The presence of certain variants may inform clinical management.

ARUP's Myeloid Malignancies Mutation and Copy Number Variation Panel by Next Generation Sequencing ([3016621](#)) uses massively parallel sequencing (also known as next generation sequencing [NGS]) to detect molecular changes (including single nucleotide variants [SNVs] and small insertions and deletions [indels]), copy number variants (CNVs) for the targeted genes, and terminal copy number-neutral loss of heterozygosity (CN-LOH). The Myeloid Malignancies Mutation Panel by Next Generation Sequencing ([2011117](#)) uses massively parallel sequencing to detect molecular changes including SNVs and small indels but does not detect CNVs or CN-LOH. Because these panels overlap, they should not be concurrently ordered. If both panels are ordered on the same specimen, test 2011117 will be canceled.

These panels are more cost-effective than using multiple single gene tests for the detection of somatic variants in myeloid malignancies and can be used to complement the morphologic and cytogenetic workup of myeloid malignancies.

## Disease Overview

### Diagnostic, Prognostic, and Treatment Issues

- Targets in this panel are relevant across the spectrum of myeloid malignancies.
- Identification of one or more clonal genetic abnormalities, variants, or patterns of variants may aid in establishing the diagnosis and classification, determining prognosis, and informing clinical management of myeloid malignancies.

## Genetics

### Genes Tested

This panel tests *ANKRD26*, *ASXL1*, *ASXL2*, *BCOR*, *BCORL1*, *BRAF*, *CALR*, *CBL*, *CBLB*, *CEBPA*, *CSF3R*, *CUX1*, *DDX41*, *DNMT1*, *DNMT3A*, *ELANE*, *ETNK1*, *ETV6*, *EZH2*, *FBXW7*, *FLT3*, *GATA1*, *GATA2*, *GNAS*, *HNRNPK*, *IDH1*, *IDH2*, *IL7R*, *JAK1*, *JAK2*, *JAK3*, *KDM6A*, *KIT*, *KMT2A*, *KRAS*, *LUC7L2*, *MPL*, *NOTCH1*, *NPM1*, *NRAS*, *NSD1*, *PHF6*, *PIGA*, *PPM1D*, *PRPF40B*, *PRPF8*, *PTPN11*, *RAD21*, *RUNX1*, *SAMD9*, *SAMD9L*, *SETBP1*, *SF3B1*, *SH2B3*, *SMC1A*, *SMC3*, *SRSF2*, *STAG2*, *STAT3*, *STAT5B*, *SUZ12*, *TET2*, *TP53*, *U2AF1*, *U2AF2*, *UBA1*, *WT1*, and *ZRSR2*.

For some genes, one or more exons of the preferred transcript are not covered by sequencing for the indicated gene. Refer to the [Genes Tested](#) table below for a full list of targeted regions and exclusions.

## Featured ARUP Testing

[Myeloid Malignancies Mutation and Copy Number Variation Panel by Next Generation Sequencing 3016621](#)

**Method:** Massively Parallel Sequencing

[Myeloid Malignancies Mutation Panel by Next Generation Sequencing 2011117](#)

**Method:** Massively Parallel Sequencing

For more information on ARUP's AML panels, which include some of the genes in this panel that are specific to AML, refer to the [Acute Myeloid Leukemia Mutation Panel by Next Generation Sequencing](#) and [Rapid Acute Myeloid Leukemia Targeted Therapy Mutation Panel Test](#) Fact Sheets.

For detection of sequence variants in the *TP53* gene alone, refer to the [Somatic \*TP53\* Mutations in Whole Blood and Bone Marrow 3017691](#) test.

For more information on ARUP's genomic microarray test offerings in oncology, refer to the [Cytogenomic Microarray, Oncology](#) Test Fact Sheet.

# Test Interpretation

## Results

Reported variants are classified into two categories:

- Tier 1: Molecular mutations, CNVs, and CN-LOH with known clinical significance in hematologic malignancies
- Tier 2: Variants of unknown clinical significance in hematologic malignancies

Clinical significance in hematologic malignancies will be described, if known.

## Reported Variants

	<a href="#">Myeloid Malignancies Mutation and Copy Number Variation Panel by Next Generation Sequencing (3016621)</a>	<a href="#">Myeloid Malignancies Mutation Panel by Next Generation Sequencing (2011117)</a>
<b>Reported</b>	Sequence variants in the preferred transcript CNVs (gains or losses) in the targeted genes Likely acquired terminal CN-LOH CNVs $\geq 5$ Mb in any gene Losses in <i>TBL1XR1</i> , <i>CD200</i> , <i>IKZF1</i> , <i>CDKN2A</i> , <i>ASMTL</i> , <i>ERG</i> , <i>ARID2</i> , and/or <i>ATM</i> Gains in <i>MYC</i> Losses between <i>FIP1L1</i> and <i>PDGFRA</i> that result in a potential fusion Any CN-LOH involving <i>TP53</i> , <i>JAK2</i> , and/or <i>CBL</i>	Sequence variants in the preferred transcript
<b>Not reported</b>	Benign or likely benign variants Likely germline or interstitial CN-LOH Due to the complexity of analysis, CNVs may not be reported in instances of stem cell transplants that present with mixed chimerism, increased genomic complexity ( $>4$ CNVs), and complex aneuploidies (e.g., hyper- or hypodiploidy) VAF for CNVs with copy numbers $>3$	Benign or likely benign variants CNVs CN-LOH

Mb, megabases; VAF, variant allele fraction

## Limitations

- Some variants may not be identified due to technical limitations in the presence of pseudogenes or in repetitive or homologous regions.
- This panel is not intended to detect minimal residual disease (MRD).
- Interpretation of panel results may be impacted if the patient has had an undisclosed allogeneic bone marrow or stem cell transplant.
- This panel does not distinguish between somatic and germline variants.
- The Myeloid Malignancies Mutation and Copy Number Variation Panel by Next Generation Sequencing ([3016621](#)) does not replace conventional cytogenetic studies or genomic microarray in the workup of hematologic malignancies.
- Neither panel detects the following types of variants:
  - Variants in regions that are not included in the preferred transcript for the targeted genes; refer to the [Genes Tested](#) table for a full list of targeted regions and exclusions
  - RNA variants
  - Gene fusions, balanced translocations, and other structural variants

## Limit of Detection

- SNVs, indels, complex variants, and *FLT3*-internal tandem duplications (*FLT3*-ITDs) are expected to be detected down to 5% VAF; variants with sufficient quality and read depth may be reported below 5% VAF.
- CNVs (gains and losses):  $>2$  Mb in approximately 30% of the sample
- CN-LOH:  $>10$  Mb in approximately 30% of the sample
  - Some areas of the genome may have a reduced sensitivity for CNVs and CN-LOH at limit of detection (LOD).

## Analytic Sensitivity

Variant Class	Analytic Sensitivity (PPA) <sup>a</sup> Estimate (%)	Analytic Sensitivity (PPA) 95% Credibility Region <sup>a</sup> (%)
SNVs	>99	99.4-100.0
Indels, duplications, complex variants	98.5	96.3-99.5
<i>FLT3</i> -ITDs	96.6	89.4-99.3
Copy number gains (>2 Mb)	91.8	86.7-95.3
Copy number losses (>2 Mb)	92.3	87.7-95.5
CN-LOH (>10 Mb)	98.1	91.5-99.8

<sup>a</sup>Genes included on this panel are a subset of a larger methods-based validation from which the PPA values are derived.

bp, base pairs; PPA, positive percent agreement

## Genes Tested

Gene	Preferred Transcript <sup>a</sup>	Excluded Exons <sup>b</sup>
<i>ANKRD26</i>	NM_014915	—
<i>ASXL1</i>	NM_015338	—
<i>ASXL2</i>	NM_018263	—
<i>BCOR</i>	NM_001123385	—
<i>BCORL1</i>	NM_021946	—
<i>BRAF</i>	NM_004333	—
<i>CALR</i>	NM_004343	—
<i>CBL</i>	NM_005188	—
<i>CBLB</i>	NM_170662	—
<i>CEBPA</i>	NM_004364	—
<i>CSF3R</i>	NM_156039	—
<i>CUX1</i>	NM_181552	24
<i>DDX41</i>	NM_016222	—
<i>DNMT1</i>	NM_001130823	5
<i>DNMT3A</i>	NM_175629	—
<i>ELANE</i>	NM_001972	—
<i>ETNK1</i>	NM_018638	—
<i>ETV6</i>	NM_001987	—
<i>EZH2</i>	NM_004456	—
<i>FBXW7</i>	NM_033632	—

Gene	Preferred Transcript <sup>a</sup>	Excluded Exons <sup>b</sup>
<i>FLT3</i>	NM_004119	—
<i>GATA1</i>	NM_002049	—
<i>GATA2</i>	NM_032638	—
<i>GNAS</i>	NM_000516	—
<i>HNRNPK</i>	NM_002140	—
<i>IDH1</i>	NM_005896	—
<i>IDH2</i>	NM_002168	—
<i>IL7R</i>	NM_002185	—
<i>JAK1</i>	NM_002227	—
<i>JAK2</i>	NM_004972	—
<i>JAK3</i>	NM_000215	—
<i>KDM6A</i>	NM_001291415	13
<i>KIT</i>	NM_000222	—
<i>KMT2A</i>	NM_001197104	—
<i>KRAS</i>	NM_004985	—
<i>LUC7L2</i>	NM_016019	—
<i>MPL</i>	NM_005373	—
<i>NOTCH1</i>	NM_017617	—
<i>NPM1</i>	NM_002520	1
<i>NRAS</i>	NM_002524	—
<i>NSD1</i>	NM_022455	—
<i>PHF6</i>	NM_001015877	—
<i>PIGA</i>	NM_002641	—
<i>PPM1D</i>	NM_003620	—
<i>PRPF8</i>	NM_006445	—
<i>PRPF40B</i>	NM_001031698	—
<i>PTPN11</i>	NM_002834	—
<i>RAD21</i>	NM_006265	—
<i>RUNX1</i>	NM_001754	—
<i>SAMD9</i>	NM_017654	—
<i>SAMD9L</i>	NM_152703	—
<i>SETBP1</i>	NM_015559	—
<i>SF3B1</i>	NM_012433	—
<i>SH2B3</i>	NM_005475	—

Gene	Preferred Transcript <sup>a</sup>	Excluded Exons <sup>b</sup>
<i>SMC1A</i>	NM_006306	—
<i>SMC3</i>	NM_005445	—
<i>SRSF2</i>	NM_003016	—
<i>STAG2</i>	NM_001042749	—
<i>STAT3</i>	NM_139276	—
<i>STAT5B</i>	NM_012448	6-9
<i>SUZ12</i>	NM_015355	1-9
<i>TET2</i>	NM_001127208	—
<i>TP53</i> <sup>c</sup>	NM_000546	—
<i>U2AF1</i>	NM_006758	—
<i>U2AF2</i>	NM_007279	—
<i>UBA1</i>	NM_003334	—
<i>WT1</i>	NM_024426	—
<i>ZRSR2</i>	NM_005089	—

<sup>a</sup>This is the transcript number used for analyzing and reporting variants. The transcript version number may change periodically and thus is not listed here. The transcript with the version number will be included on the patient's report if a variant is detected in the gene.

<sup>b</sup>Noncoding exons are not analyzed, except for regions containing known clinically relevant variants in the *ANKRD26* 5'UTR and *NOTCH1* 3'UTR. In addition, coding exons noted here are not sequenced due to technical limitations of the assay.

<sup>c</sup>For detection of sequence variants in the *TP53* gene alone, refer to the [Somatic \*TP53\* Mutations in Whole Blood and Bone Marrow 3017691](#) test.

ARUP Laboratories is a nonprofit enterprise of the University of Utah and its Department of Pathology, 500 Chipeta Way, Salt Lake City, UT 84108  
(800) 522-2787 | (801) 583-2787 | aruplab.com | arupconsult.com