

Rett Syndrome and *MECP2*-Related Disorders

Indications for Ordering

- Confirm clinical diagnosis of Rett syndrome or *MECP2*-related disorder
- Determine cause of severe neonatal encephalopathy or intellectual disability in males
- Rule out *MECP2* variant in individuals with clinical features of Angelman syndrome who lack a molecular abnormality involving 15q11.2-13

Test Description

- Polymerase chain reaction followed by bidirectional sequencing of the *MECP2* coding region (exons 1-4) and intron/exon boundaries
- Multiplex ligation-dependent probe amplification to detect large deletions/duplications in *MECP2* coding regions (exons 1-4)

Tests to Consider

Testing strategy

Consultation with a genetic counselor is recommended to plan the optimal *MECP2* genetic testing sequence

- Sequencing
 - Consider for **females** with
 - Classic or atypical Rett syndrome
 - Nonspecific intellectual disability
 - Autism
 - Clinically suspected, but molecularly unconfirmed, Angelman syndrome
 - Consider for **males** with
 - Unexplained neonatal encephalopathy
 - Nonspecific intellectual disability
 - Phenotype of Rett syndrome (classic or atypical)
 - Manic-depressive psychosis, pyramidal signs, Parkinsonian and macroorchidism (PPM-X syndrome)
 - Clinically suspected, but molecularly unconfirmed, Angelman syndrome

- Deletion/duplication
 - Consider for **females** when suspicion for classic or atypical Rett syndrome remains despite negative *MECP2* sequencing result
 - Consider for **males** with suspected *MECP2* duplication syndrome
 - Severe developmental delay/intellectual disability and a history of recurrent infections, swallowing dysfunction, hypotonia, spasticity, limited speech, and limited ambulation

Primary tests

[Rett Syndrome \(*MECP2*\), Sequencing and Deletion/Duplication 0051614](#)

- Preferred initial test for females with a phenotype of classic or atypical Rett syndrome
- Detects up to 95% of *MECP2* pathogenic variants

[Rett Syndrome \(*MECP2*\), Full Gene Sequencing 0051378](#)

- Acceptable first-tier diagnostic test for females with a phenotype of classic or atypical Rett syndrome

Related tests

[Familial Mutation, Targeted Sequencing 2001961](#)

- Useful when a pathogenic familial variant identifiable by sequencing is known

[CDKL5-Related Disorders \(*CDKL5*\) Sequencing and Deletion/Duplication 2004935](#)

- May be useful for individuals with early-onset seizures and intellectual disability when *MECP2* full-gene analysis fails to detect a pathogenic variant

Disease Overview

Rett syndrome

Incidence

- 1/10,000-15,000 female births

Clinical diagnostic criteria (Neul, 2010)

- Required for classic (typical) Rett syndrome
 - Period of regression followed by recovery or stabilization
 - All main criteria
 - Supportive criteria are not required, although are often present in classic Rett
- Required for atypical Rett
 - Period of regression followed by recovery or stabilization
 - At least 2 of the 4 main criteria
 - 5 of 11 supportive criteria
- Main criteria
 - Partial or complete loss of acquired purposeful hand skills
 - Partial or complete loss of acquired spoken language
 - Gait abnormalities (impaired or absence of ability)
 - Stereotypic hand movements (wringing/squeezing, clapping/tapping, washing/rubbing, and mouthing of hands)
- Supportive criteria
 - Breathing disturbances when awake
 - Bruxism when awake
 - Impaired sleep pattern
 - Abnormal muscle tone
 - Peripheral vasomotor disturbances
 - Scoliosis/kyphosis
 - Growth retardation
 - Small, cold hands and feet
 - Inappropriate laughing/screaming spells
 - Diminished response to pain
 - Intense eye communication

Clinical phenotype – highly variable

Influenced by the patient's sex, age, the specific *MECP2* variant and, in females, the X-inactivation pattern

MECP2 duplication syndrome

Incidence

- 1-2% of males with severe neonatal encephalopathy or moderate to severe intellectual disability
- Female carriers tend to be asymptomatic

Clinical phenotype – suspected in males with

- Severe intellectual disability
- Limited/lack of speech
- Early-onset hypotonia
- Progressive spasticity
- Seizures
- Recurrent respiratory infections
- Limited/lack of ambulation

Genetics

Gene – *MECP2*

Inheritance – X-linked dominant

- 99.5% of cases are sporadic

Penetrance

- Nearly 100% in females
- Female variant carriers may be asymptomatic as a result of highly skewed X-inactivation
- Disruption of *MECP2* in males is usually lethal

Variants – >350 are known

- Sequence variants
 - Classic Rett syndrome – detected in ~80% of affected individuals
 - Atypical Rett syndrome – detected in ~40% of affected individuals
 - Majority of *MECP2* sequence variants occur in exon 4
- Deletions/duplications
 - Large deletions most commonly involve exon 3 and/or exon 4 of *MECP2*
 - Large deletions account for ~8% of classic Rett syndrome and ~3% of atypical Rett syndrome
 - Males with whole-gene duplications are clinically affected

Test Interpretation

Sensitivity/specificity

- Clinical sensitivity
 - Classic Rett syndrome
 - Sequencing – ~80% (Fukuda, 2005; Li, 2007)
 - Deletions/duplications – ~8-10% (Archer, 2006; Li, 2007)
 - Atypical Rett
 - Sequencing – ~40% (Fukada, 2005; Li, 2007)
 - Deletions/duplications – ~3% (Archer, 2006)
- Analytical sensitivity
 - Sequencing – 99%
 - Deletions/duplications – 90%
- Analytical specificity
 - Sequencing – 99%
 - Deletions/duplications – 98%

Results

- Positive – pathogenic variant identified
 - Diagnosis confirmed
- Negative – no pathogenic variant identified
 - Possibility of an *MECP2*-related disorder is decreased, but not excluded
- Inconclusive – gene variant detected
 - Unknown whether variant is pathogenic or benign

Limitations

- Diagnostic errors can occur due to rare sequence variations
- Not detected
 - Breakpoints of large deletions/duplications
 - Regulatory region and deep intronic variants

References

- Archer HL, Whatley SD, et al. Gross rearrangements of the MECP2 gene are found in both classical and atypical Rett syndrome patients. *J Med Genet.* 2006;43(5):451-456
- Fukuda T, Yamashita Y, et al. Methyl-CpG binding protein 2 gene (MECP2) variations in Japanese patients with Rett syndrome: pathological mutations and polymorphisms. *Brain Dev.* 2005;27(3):211-217
- Li MR, Pan H, et al. MECP2 and CDKL5 gene mutation analysis in Chinese patients with Rett syndrome. *J Hum Genet.* 2007;52(1):38-47
- Neul JL, Kaufmann WE, et al. Rett syndrome: revised diagnostic criteria and nomenclature. *Ann Neurol.* 2010;68(6):944-950