

TEST CHANGE

Gamma Glutamyl Transferase, Serum or Plasma

0020009, GT

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube or plasma separator tube. Also acceptable: **Green (lithium heparin)** ~~Pink (K₂EDTA)~~.

Specimen Preparation: Allow serum tube to clot completely at room temperature. Separate serum or plasma from cells within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.2 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Grossly hemolyzed specimens.

Remarks:

Stability: After separation from cells: Ambient: 1 week; Refrigerated: 1 week; Frozen: 2 months

Methodology: Quantitative Enzymatic Assay

Performed: Sun-Sat

Reported: Within 24 hours

Note:

CPT Codes: 82977

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report (reports may vary based on instrumentation)