

HOTLINE: Effective February 22, 2022

0020224

Myoglobin, Serum

MYOG-S

Specimen Required: Collect: Plain red or serum separator tube. Also acceptable: Green (lithium heparin), lavender (K₃EDTA or K₂EDTA) or pink (K₂EDTA).

Specimen Preparation: Allow specimen to clot completely at room temperature. Separate serum or plasma from cells ASAP or within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.2 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Grossly hemolyzed specimens.

Stability (collection to initiation of testing): After separation from cells: Ambient: 8 days; Refrigerated: 2 weeks; Frozen: 1 year

Reference Interval:

Effective February 22, 2022

Male: less than or equal to 72 ng/mL

Female: less than or equal to 58 ng/mL

HOTLINE NOTE: There is a numeric map change associated with this test.

Change the numeric map for component 0020224, Myoglobin Serum from XXXXXX to XXXXXXXX.