

## **TEST CHANGE**

Lactate Dehydrogenase, Isoenzymes 0020413, LDISO	
Specimen Requirements:	
Patient Preparation:	
Collect:	Serum <u>separator tube</u> Separator Tube (SST) or <u>plain red</u> Plain <del>Red</del> .
Specimen Preparation:	Do not refrigerate or freeze. Allow specimen to clot completely at room temperature. Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP <u>standard</u> <u>transport tube.Standard Transport Tube.</u> (Min: 0.6 mL)
Transport Temperature:	Room temperature.
Unacceptable Conditions:	Hemolyzed specimens.
Remarks:	
Stability:	After separation of cells: Ambient: 1 week; Refrigerated: Unacceptable; Frozen: Unacceptable
Methodology:	Quantitative Enzymatic Assay/Electrophoresis
Performed:	Sun-Sat
Reported:	1- <u>3</u> 2 days
Note:	LD-1 and LD-2 are elevated in hemolyzed specimens and serum which has not been separated from cells. LD-3, LD-4, and LD-5 are labile at low temperatures, and are erroneously low in specimens that have been refrigerated or frozen.
CPT Codes:	83625; 83615
New York DOH Approval Status:	This test is New York DOH approved.
Interpretive Data:	
Reference Interval:	
LD-1: 14-27% LD-2: 29-42% LD-3: 18-30% LD-4: 8-15% LD-5: 6-23%	



Lactate Dehydrogenase, Total: 0-1 month: 200-465 U/L 2-17 months: 200-450 U/L 18 months-10 years: 165-430 U/L 11-16 years: 127-287 U/L 17 years and older: 105-230 U/L