

TEST CHANGE

Procalcitonin

0020763, PCT

Specimen Requirements:

Patient Preparation: The same specimen type (serum, plasma) should be used throughout the patient's clinical course.

Collect: Plasma ~~separator tube~~ Separator Tube (PST) or serum separator tube ~~Serum Separator Tube~~ (SST).

Specimen Preparation: ~~For~~ Allow serum ~~specimens, ensure that complete~~ to sit for 15-20 minutes for proper clot formation has taken place prior to centrifugation. If and to ensure the specimen is centrifuged before complete clot formation, the presence ~~absence~~ of fibrin may cause erroneous results. The use of plasma is recommended for rapid turnaround of results. For accurate results, in the ~~serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter, which can interfere with this assay.~~ Separate from cells ASAP or within 2 hours of collection. Transfer 2 mL serum or plasma to an ARUP standard transport tube. ~~Standard Transport Tube.~~ (Min: 0.3 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Specimens collected in citrate anticoagulant. Specimens that are heat-inactivated, pooled, grossly hemolyzed, contain obvious microbial contamination or fungal growth should not be used.

Remarks:

Stability: After separation from cells: Ambient: 24 hours; Refrigerated: 5 days; Frozen: 15 days

Methodology: Quantitative Chemiluminescent Immunoassay (CLIA)

Performed: Sun-Sat

Reported: Within 24 hours

Note: Procalcitonin levels below 0.50 ng/mL do not exclude an infection, because localized infections (without systemic signs) may also be associated with such low levels.

CPT Codes: 84145

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

~~Procalcitonin > A correction has been applied to optimize cutoffs established for the BRAHMS PCT sensitive KRYPTOR assay.~~

~~Procalcitonin > 2.00 ng/mL: Procalcitonin levels above 2.00 ng/mL on the first day of ICU admission represent a high risk for progression to severe sepsis and/or septic shock.~~

~~Procalcitonin < 0.50 ng/mL: Procalcitonin levels below 0.50 ng/mL on the first day of ICU admission represent a low risk for progression to severe sepsis and/or septic shock.~~

~~If the procalcitonin measurement is performed shortly after the systemic infection process has started (usually less than 6 hours), these values may still be low. As various ~~noninfectious~~ ~~non-~~ ~~infectious~~ conditions are known to induce procalcitonin as well, procalcitonin levels between 0.50 ng/mL and 2.00 ng/mL should be reviewed carefully to take into account the specific clinical background and condition(s) of the individual patient.~~

Reference Interval:

Less than 0.07 ng/mL