

TEST CHANGE

APC Resistance Profile

0030127, APC RST

Specimen Requirements:

Patient Preparation:

Collect: Lt. blue (sodium citrate). ~~Special~~Refer to Specimen Collection and Handling Hemostasis/Thrombosis Specimens guide located at <https://www.aruplab.com/Specimen-Handling/SpecialSpecimenCollection/Hemostasis-Thrombosis.pdf> for hemostasis/thrombosis specimen handling guidelines.

Specimen Preparation: Transfer 1.5 mL platelet-poor plasma to an ARUP standard transport tube. ~~Standard Transport Tube~~ (Min: 1 mL)

Transport Temperature: CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.

Unacceptable Conditions: Serum. EDTA plasma, clotted or hemolyzed specimens.

Remarks:

Stability: Ambient: 4 hours; Refrigerated: Unacceptable; Frozen ~~at -20 Degrees C~~: 3 months; ~~Frozen at -70 Degrees C~~: 6 months

Methodology: Electromagnetic Mechanical Clot Detection

Performed: Mon-Sat

Reported: 1-4 days

Note:

CPT Codes: 85307

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Ratios less than 2.00 suggest APC resistance. This method uses factor V deficient plasma; therefore, APC resistance due to a nonfactor V mutation will not be detected. Extreme factor V deficiency or presence of direct oral anticoagulants (DOACs) may cause an unreliable ratio.

Reference Interval:

2.00 or greater

