
0030143 Heparin Anti-Xa, Unfractionated**HA UN**

Performed: Sun-Sat
Reported: 1-2 days

Specimen Required: Collect: Lt. blue (sodium citrate). Refer to Specimen Handling at aruplab.com for hemostasis/thrombosis specimen handling guidelines.
Specimen Preparation: Centrifuge specimen within one hour of collection. Transport 2 mL platelet-poor plasma. (Min: 1 mL)
Storage/Transport Temperature: Frozen.
Remarks: This test cannot be used to quantitate anticoagulants other than unfractionated heparin. This includes, but is not limited to, direct oral anticoagulants and Fondaparinux (Arixtra).
Unacceptable Conditions: Serum. EDTA, oxalate, heparin, or plasma separator tubes. Hemolyzed or clotted specimens.
Stability (collection to initiation of testing): After separation from cells: Ambient: **Unacceptable**; Refrigerated: **Unacceptable**; Frozen: 1 month
University of Utah Clients: After separation from cells: Ambient: 4 hours; Refrigerated: 8 hours; Frozen: 2 weeks