

0030192

APC Resistance Profile with Reflex to Factor V Leiden**APC R**

Specimen Required: Collect: Light Blue (Sodium Citrate) **AND** Lavender (EDTA), Pink (K₂EDTA), or Yellow (ACD Solution A or B). Refer to Specimen Handling at aruplab.com for hemostasis/thrombosis specimen handling guidelines.

Specimen Preparation: Transport 1.5 mL platelet-poor plasma **AND** 3 mL whole blood. (Min: 1 mL/each)

Storage/Transport Temperature: **Plasma:** CRITICAL FROZEN. **Separate specimens must be submitted when multiple tests are ordered.**

Whole Blood: Refrigerated.

Unacceptable Conditions: Serum, clotted or hemolyzed specimens. **Frozen specimens in glass collection tubes.**

Stability (collection to initiation of testing): **Plasma:** Ambient: 4 hours; Refrigerated: Unacceptable; Frozen at -20°C: 3 months; Frozen at -70°C: 6 months

Whole Blood: Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable