

HOTLINE: Effective February 22, 2022

0050149 Complement Components 3 and 4

C3C4

Specimen Required: Collect: Serum separator tube, Plasma separator tube, Green (Lithium heparin)

Specimen Preparation: Allow specimen to clot at room temperature. Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube and freeze. (Min: 0.6 mL)

 $\underline{Storage/Transport\ Temperature:}\ Frozen.$

<u>Unacceptable Conditions:</u> Specimens left to clot at refrigerated temperature. Specimens exposed to repeated freeze/thaw cycles. Stability (collection to initiation of testing): Ambient: 48 hours; Refrigerated: 8 days; Frozen: 2 weeks (avoid repeated freeze/thaw

Reference Interval:

Test Number	Components	Reference Interval	
0050150	Complement Component 3	Effective February 22, 2022	
		0-30 days: 59-121 mg/dL 1 month: 55-129 mg/dL 2 months: 61-155 mg/dL 3 months: 67-136 mg/dL 4 months: 65-182 mg/dL 5 months: 67-174 mg/dL 6 months: 77-179 mg/dL	7-8 months: 78-173 mg/dL 9-11 months: 76-187 mg/dL 1 year: 87-181 mg/dL 2 years: 84-177 mg/dL 3-4 years: 80-178 mg/dL 5-11 years: 80-160 mg/dL 12-17 years: 82-163 mg/dL 18 years and older: 90-180 mg/dL
0050155	Complement Component 4		
0000133	Compensar Component	0-30 days: 8-30 mg/dL 1 month: 9-33 mg/dL 2 months: 9-37 mg/dL 3 months: 10-35 mg/dL 4 months: 10-49 mg/dL 5 months: 9-48 mg/dL 6 months: 12-55 mg/dL	7-8 months: 13-48 mg/dL 9-11 months: 16-51 mg/dL 1 year: 16-52 mg/dL 2-4 years: 12-47 mg/dL 5-11 years: 13-44 mg/dL 12-17 years: 14-41 mg/dL 18 years and older: 10-40 mg/dL