

TEST CHANGE

Varicella-Zoster Virus Antibody, IgG

0050167, VZE

Specimen Requirements:

Patient Preparation:

Collect: Serum ~~separator tube~~ **Separator Tube** (SST).

Specimen Preparation: Allow specimen to clot completely at room temperature. Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum an ARUP ~~standard transport tube~~ **Standard Transport Tube**. (Min: 0.5 mL) Parallel testing is preferred and convalescent specimens must be received within 30 days from receipt of the acute specimens.

Transport Temperature: Refrigerated.

Unacceptable Conditions: Contaminated, heat-inactivated, grossly hemolyzed, lipemic, or severely icteric specimens.

Remarks: Label specimens plainly as acute or convalescent.

Stability: After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year (~~a~~**A**void repeated freeze/thaw cycles)

Methodology: Semi-Quantitative Chemiluminescent Immunoassay

Performed: Sun-Sat

Reported: Within 24 hours

Note: For CSF specimens, refer to Varicella-Zoster Virus Antibody, IgG, CSF (ARUP test code 0054444).

CPT Codes: 86787

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

The best evidence for current infection is a significant change on two appropriately timed specimens, where both tests are done in the same laboratory at the same time.

Reference Interval:

Effective **October 7, 2024**

~~Q~~ **February 18, 2020**

~~134.9 S/COIV~~ or less: Negative - No significant level of detectable varicella-zoster IgG antibody.

~~135.0-164.9 IV: Equivocal - Repeat testing in 10-14 days may be helpful.~~

~~165.0 S/COIV~~ or greater: Positive - IgG antibody to varicella-zoster detected, which may indicate a current or past varicella-zoster infection.

HOTLINE NOTE: There is a numeric map change associated with this test. Refer to the Hotline Test Mix for interface build information.

HOTLINE NOTE: There is a unit of measure change associated with this test. Refer to the Hotline Test Mix for interface build information.