

## TEST CHANGE

### Haptoglobin

0050280, HAPTO

#### Specimen Requirements:

Patient Preparation: Fasting specimen preferred.

Collect: Plasma separator tube or serum separator tube. Also acceptable: Green (lithium heparin), or pink (K2EDTA).

Specimen Preparation: Allow specimen to clot completely at room temperature. Separate serum or plasma from cells ASAP or within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP standard transport tube. ~~Standard Transport Tube~~. (Min: 0.5 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Grossly hemolyzed specimens.

#### Remarks:

Stability: After separation from cells: Ambient: 3 months; Refrigerated: 8 months; Frozen: 3+ months

Methodology: Quantitative Immunoturbidimetry

Performed: Sun-Sat

Reported: Within 24 hours

#### Note:

CPT Codes: 83010

New York DOH Approval Status: This test is New York DOH approved.

#### Interpretive Data:

#### Reference Interval:

30-200 mg/dL