

TEST CHANGE

Immunoglobulin E

0050345, IGE

Specimen Requirements:

Patient Preparation:

Collect:

Serum separator tube (SST) or plasma separator tube (PST).
Also acceptable: Green (sodium or lithium heparin), lavender
(K2EDTA), or pink (K2EDTA).

Specimen Preparation:

Separate from cells ASAP or within 2 hours of collection.
Transfer 1 mL serum or plasma to an ARUP standard transport
tube. (Min: 0.3 mL)

Transport Temperature:

Refrigerated.

Unacceptable Conditions:

Hemolyzed, icteric, or lipemic specimens. Postmortem
samples.

Remarks:

Stability:

After separation from cells: Ambient: 48 hours; Refrigerated: 2
weeks; Frozen: 1 year

Methodology:

Quantitative ImmunoCAP Fluorescent Enzyme Immunoassay

Note:

CPT Codes:

82785

New York DOH Approval Status:

This test is New York DOH approved.

Interpretive Data:

Age	Reference Interval
0-5 months	13 kU/L or less
6-12 months	34 kU/L or less
1-2 years	97 kU/L or less
3 years	199 kU/L or less
4-6 years	307 kU/L or less
7-8 years	403 kU/L or less
9-12 years	696 kU/L or less
13-15 years	629 kU/L or less
16-17 years	537 kU/L or less
18 years and older	214 kU/L or less

Inserted Cells

Inserted Cells

Reference Interval:

Effective November 17, 2014

Test Number	Components	Reference Interval
	Age	Reference Interval
	0-5 months	13 kU/L or less
	6-12 months	34 kU/L or less
	1-2 years	97 kU/L or less
	3 years	199 kU/L or less
	4-6 years	307 kU/L or less
	7-8 years	403 kU/L or less
	9-12 years	696 kU/L or less
	13-15 years	629 kU/L or less
	16-17 years	537 kU/L or less
	18 years and older	214 kU/L or less
	Immunoglobulin E	

Age	Reference Interval (kU/L)
0-5 months	13 or less
6-12 months	34 or less
1-2 years	97 or less
3 years	199 or less
4-6 years	307 or less
7-8 years	403 or less
9-12 years	696 or less
13-15 years	629 or less
16-17 years	537 or less

Deleted Cells

Inserted Cells

