

HOTLINE: Effective February 16, 2021

0050576	Immunoglobulin G Subclass 4	IGG4
----------------	------------------------------------	-------------

Methodology: Quantitative **Immunoturbidimetry**
Performed: Sun-Sat
Reported: 1-3 days

Specimen Required: Collect: Serum separator tube.

Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: **0.5** mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: **Grossly hemolyzed or lipemic specimens**

Stability (collection to initiation of testing): After separation from cells: Ambient: **Unacceptable**; Refrigerated: **14** days; Frozen: 6 months

Reference Interval:

Effective February 16, 2021

Age	Reference Interval
0-2 years	1-34 mg/dL
3-4 years	1-65 mg/dL
5-9 years	0-168 mg/dL
10-14 year	1-103 mg/dL
15-18 years	2-170 mg/dL
19 years and older	1-123 mg/dL