

HOTLINE: Effective **November 15, 2021**

0051627

Epstein-Barr Virus Antibody to Viral Capsid Antigen, IgG and IgA

EBV PAN 3

Specimen Required: Collect: Serum Separator Tube (SST).

Specimen Preparation: Allow specimen to clot completely at room temperature. Separate from cells ASAP or within 2 hours of collection. Transfer 2 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL) Parallel testing is preferred and convalescent specimens **must** be received within 30 days from receipt of the acute specimens.

Storage/Transport Temperature: Refrigerated.

Remarks: **Label specimens plainly as acute or convalescent.**

Unacceptable Conditions: Contaminated, heat-inactivated, or grossly hemolyzed specimens.

Stability (collection to initiation of testing): After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 **month** (Avoid repeated freeze/thaw cycles).

Reference Interval:

Test Number	Components	Reference Interval	
0050235	Epstein-Barr Virus Antibody to Viral Capsid Antigen, IgG	Effective February 19, 2013	
		17.9 U/mL or less	Not Detected
		18.0-21.9 U/mL	Indeterminate - Repeat testing in 10-14 days may be helpful.
		22.0 U/mL or greater	Detected
0051626	Epstein-Barr Virus Antibody to Viral Capsid Antigen, IgA	Effective November 15, 2021	
		8 U or less	Not Detected
		9-11 U	Indeterminate - Repeat testing in 10-14 days may be helpful.
		12 U or greater	Detected

HOTLINE NOTE: There is a numeric map change associated with this test.

Change the numeric map for component 0051626, EBV Antibody To Viral Capsid Antigen IgA from XX.X to **XX**.

There is a unit of measure change associated with this test.

Change the unit of measure for component 0051626, EBV Antibody To Viral Capsid Antigen IgA from U/L to **U**.