
0055593

Hepatitis C Virus Genotype by Sequencing

HEPCGENO

Specimen Required: Collect: Lavender (EDTA), Pink (K₂EDTA), Plasma Preparation Tube (PPT) or Serum Separator Tube (SST).
Specimen Preparation: Separate from cells within 6 hours of collection. Transfer 2 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.5 mL)
Storage/Transport Temperature: Frozen.
Remarks: Please submit most recent viral load and test date, if available.
Unacceptable Conditions: Heparinized specimens.
Stability (collection to initiation of testing): After separation from cells: Ambient: Unacceptable; Refrigerated: 72 hours; Frozen: 4 months