

TEST CHANGE

Citric Acid, Serum or Plasma

0080011, CITRIC S

Specimen Requirements:

Patient Preparation:

Collect: Plain red or serum separator tube (SST). Also acceptable: Green (sodium heparin), lavender (EDTA), or plasma preparation tube (PPT) ~~white (potassium EDTA)~~.

Specimen Preparation: Transfer 2 mL serum or plasma to an ARUP standard transport tube ~~Standard Transport Tube~~. (Min: 1 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Refrigerated. Also acceptable: Frozen.

Unacceptable Conditions: Room temperature specimens.

Remarks:

Stability: Ambient: 8 hours; Refrigerated: 17 ~~18~~ days; Frozen: 17 days ~~1 month~~

Methodology: Quantitative Spectrophotometry/Enzymatic Assay

Performed: Varies

Reported: 3-12 ~~11~~ days

Note:

CPT Codes: 82507

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Reference Interval:

By report