

TEST CHANGE

Porphobilinogen (PBG), Urine

0080260, PBGQT

Specimen Requirements:

Patient Preparation:

Collect: Random or 24-hour urine. Refrigerate 24-hour specimens during collection.

Specimen Preparation: Protect from light. Transfer 28 mL aliquot from a random or well-mixed 24-hour collection to ARUP amber transport tubes. ~~Amber Transport Tubes~~. (Min: 13.5 mL) Record total volume and collection time interval on transport tube and test request form.

Transport Temperature: Frozen.

Unacceptable Conditions: Body fluids other than urine.

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: 1 week~~4 days~~; Frozen: 1 month

Methodology: Quantitative Ion Exchange Chromatography/Spectrophotometry

Performed: Sun-Sat~~Mon-Fri~~

Reported: 1-5~~4~~ days

Note: Appropriate test to rule out acute intermittent porphyria (AIP) and other acute attack types of porphyrias associated with neurologic and/or psychiatric symptoms.

CPT Codes: 84110

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Porphobilinogen (PBG), Urine

Results for random urine specimens are normalized to creatinine (CRT) concentration and reported as a ratio of amounts (millimoles of PBG/mole of creatinine).

Porphobilinogen (PBG) in a random urine specimen is used to evaluate an attack of acute porphyria. Slight increases in urinary PBG are associated with acute porphyrias other than acute

intermittent porphyria (AIP) and may indicate a resolving or treated acute porphyria.

Urinary PBG in excess of two times the upper reference limit is consistent with acute porphyria.

Per 24h calculations are provided to aid interpretation for collections with a duration of 24 hours and an average daily urine volume. For specimens with notable deviations in collection time or volume, ratios of analytes to a corresponding urine creatinine concentration may assist in result interpretation.

~~This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.~~

Reference Interval:

Test Number	Components	Reference Interval		
	Creatinine, Urine - per 24h			
		Age	Male (mg/d)	Female (mg/d)
		3-8 years	140-700	140-700
		9-12 years	300-1300	300-1300
		13-17 years	500-2300	400-1600
		18-50 years	1000-2500	700-1600
		51-80 years	800-2100	500-1400
		81 years and older	600-2000	400-1300
	Porphobilinogen, Urine - per 24h	0.4 - 1.50-11.0 μmol/d		
	Porphobilinogen (PBG) , Urine - <u>ratio to CRT</u> per volume	0.0 - 0.2 mmol/mol CRT-8.8 μmol/L		

HOTLINE NOTE: There is a component change associated with this test. One or more components have been added or removed. Refer to the Hotline Test Mix for interface build information.