

TEST CHANGE

Librium and Nordiazepam

0090148, LIB

Specimen Requirements:

Patient Preparation: Timing of specimen collection: ~~Pre-dose~~ **Pre-dose** (trough) draw ~~at~~ **At** steady state concentration.

Collect: Gray (~~potassium oxalate/sodium fluoride~~), ~~Potassium Oxalate/Sodium Fluoride~~. Also acceptable: Plain **red, green** (~~sodium heparin~~), ~~lavender~~ **Red, Green (Sodium Heparin)**, ~~Lavender~~ (K2 or K3EDTA) or pink (K2EDTA).

Specimen Preparation: Separate serum or plasma from cells ASAP or within 2 hours of collection. Transfer 2 mL serum or plasma to an ARUP ~~standard transport tube~~, **Standard Transport Tube**. (Min: 1 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Gel separator tubes. Plasma or whole blood collected in light blue (sodium citrate). Hemolyzed specimens.

Remarks:

Stability: After separation from cells: Ambient: 1 week; Refrigerated: 2 weeks; Frozen: 3 years (Avoid repeated freeze/thaw cycles)

Methodology: Quantitative Liquid Chromatography-Tandem Mass Spectrometry

Performed: Tue, Fri

Reported: 1-~~7~~**5** days

Note:

CPT Codes: 80346 (Alt code: G0480)

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Adverse effects may include drowsiness, ataxia, nausea, and constipation.

Reference Interval:

Effective November 16, 2015

Components	Therapeutic Range
Librium	500-3000 ng/mL - Dose (Adult): 5-100 mg Toxic: Greater than 5000 ng/mL
Nordiazepam	100-1500 ng/mL - Based on normal dosages. Toxic: Greater than 2500 ng/mL