

TEST CHANGE

Cortisol/Cortisone Urine Free by LC-MS/MS

0092100, CORTURATIO

Specimen Requirements:

Patient Preparation:

Collect: 24-hour or random urine. Refrigerate 24-hour specimen during

collection.

Specimen Preparation: Transport one 4 mL aliquot of urine. (Min: 1 mL) Record total

volume and collection time interval on transport tube and test

Effective Date: February 20, 2024

request form.

Transport Temperature: Refrigerated.

Unacceptable Conditions: Room temperature specimens. Acidified specimens or

specimens with preservatives.

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: 1 week; Frozen: 1 month

Methodology: Quantitative High Performance Liquid Chromatography-

Tandem Mass Spectrometry

Performed: Sun-Sat

Reported: 1-54 days

Note: Reference intervals for children are based on literature from

Taylor R.L. et al., Validation of a High-Throughput Liquid Chromatography-Tandem Mass Spectrometry Method for Urine Cortisol and Cortisone. Clinical Chemistry 2002; 48:1511-19. *The ratio of the concentrations of cortisol to cortisone will not

be evaluated if the cortisol concentration is less than 5

microg/L.

CPT Codes: 82530; 83789

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Per 24h calculations are provided to aid interpretation for collections with a duration of 24 hours and an average daily urine volume. For specimens with notable deviations in collection time or volume, ratios of analytes to a corresponding urine creatinine concentration may assist in result interpretation.



Effective Date: February 20, 2024

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

Reference Interval:

Test Number	Components	Reference Interval			
	Creatinine, Urine - per 24h				
		Age	Male (mg/d)	Female (mg/d)	
		3-8 years	140-700	140-700	
		9-12 years	300-1300	300-1300	
		13-17 years	500-2300	400-1600	
		18-50 years	1000-2500	700-1600	
		51-80 years	800-2100	500-1400	
		81 years and older	600-2000	400-1300	
	Cortisol/Cortisone Ratio	18 years and	d older: 0.15-0.50		
	Cortisol, Urine Free - ratio to CRT				
		Age	Male (ug/g CRT)	Female (ug/g CRT)	
		Prepubertal	Less than 25	Less than 25	
		18 years and older	Less than 32	Less than 24	
		Pregnancy	Not Applicable	Less than 59	
	Cortisol, Urine Free - per 24h				
		Age	Male (ug/24 h)	Female (ug/24 h)	
		3-8 years	Less than or equal to 18	Less than or equal to 18	
		9-12 years	Less than or equal to 37	Less than or equal to 37	
		13-17 years	Less than or equal to 56	Less than or equal to 56	
		18 years and older	Less than or equal to 60	Less than or equal to 45	