

**TEST CHANGE**

**Drug Panel 7, Urine - Screen with Reflex to Confirmation/Quantitation**

0092184, CDASU 7

**Specimen Requirements:**

**Patient Preparation:**

**Collect:** Random urine.

**Specimen Preparation:** Transfer 8 mL urine with no additives or preservative to ARUP Standard Transport Tubes. (Min: 4 mL)

**Transport Temperature:** Refrigerated.

**Unacceptable Conditions:** Specimens exposed to repeated freeze/thaw cycles.

**Remarks:**

**Stability:** Ambient: 1 week; Refrigerated: 1 month; Frozen: 1 month

**Methodology:** Qualitative Enzyme Multiplied Immunoassay Technique/  
Quantitative Gas Chromatography-Mass Spectrometry/Quantitative Liquid Chromatography-Tandem Mass Spectrometry

**Performed:** Sun-Sat

**Reported:** 1-4 days

**Note:** If the specimen screens positive, then Confirmation/Quantitation by GC/MS and/or LC-MS/MS will be added to confirm result. Additional charges apply.

**CPT Codes:** 80307; if reflexed, add 80325; 80345; 80346; 80349; 80353; 80359; 80361; 80365; 83992 (Reflexed Alt Code: G0480)

**New York DOH Approval Status:** This test is New York DOH approved.

**Interpretive Data:**

The absence of expected drug(s) and/or drug metabolite(s) may indicate non-compliance, inappropriate timing of specimen collection relative to drug administration, poor drug absorption, diluted/adulterated urine, or limitations of testing. The concentration at which the screening test can detect a drug or metabolite varies within a drug class. Specimens for which drugs or drug classes are detected by the screen are reflexed to a second, more specific technology (GC/MS and/or LC-MS/MS). The concentration value must be greater than or equal to the cutoff to be reported as positive. Interpretive questions should be directed to the laboratory.

For medical purposes only; not valid for forensic use.

Oxycodone results are reported with the opiates results. MDMA results are reported with the amphetamines results.

Reference Interval:

~~Effective August 17, 2015~~

**Drugs Covered and Cutoff Concentrations**

Drugs/Drug Classes	Screen
THC (Cannabinoids)	<del>50</del> 20 ng/mL
Cocaine	150 ng/mL
Opiates	300 ng/mL
Oxycodone	100 ng/mL
Phencyclidine	25 ng/mL
Amphetamines	300 ng/mL
MDMA (Ecstasy)	500 ng/mL
Barbiturates	200 ng/mL
Benzodiazepines	200 ng/mL