

TEST CHANGE

Arsenic, Whole Blood

0099045, ARS B

Specimen Requirements:

Patient Preparation: Diet, medication, and nutritional supplements may introduce interfering substances. Patient should be encouraged to discontinue nutritional supplements, vitamins, minerals, and non-essential over-the-counter medications (upon the advice of their physician) and avoid shellfish and seafood for 48 to 72 hours.

Collect: Royal blue (K2EDTA) or royal blue (NaHep).

Specimen Preparation: Transport 3 or 6 mL whole blood in the original collection tube. (Min: 0.5 mL)

Transport Temperature: Room temperature. Also acceptable: Refrigerated.

Unacceptable Conditions: Specimens collected in tubes other than royal blue (K2EDTA) or royal blue (NaHep). Specimens transported in containers other than royal blue (K2EDTA) or royal blue (NaHep) tube or trace element-free transport tube. Clotted specimens.

Remarks:

Stability: Ambient: Indefinitely; Refrigerated: Indefinitely; Frozen: Unacceptable

Methodology: Quantitative Inductively Coupled Plasma-Mass Spectrometry (ICP-MS)

Performed: Sun-Sat

Reported: 1-3 days

Note:

CPT Codes: 82175

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

~~Elevated results may be due to skin or collection-related contamination, including the use of a noncertified metal-free collection/transport tube. If contamination concerns exist due to elevated levels of blood arsenic, confirmation with a second specimen collected in a certified metal-free tube is recommended.~~



Potentially toxic ranges for blood arsenic: Greater than or equal to 600 $\mu\text{g}/\text{g}/\text{L}$.

Blood arsenic is for the detection of recent exposure poisoning only. Blood arsenic levels in healthy subjects vary considerably with exposure to arsenic in the diet and the environment. A 24-hour urine arsenic is useful for the detection of chronic exposure.

Elevated results may be due to skin- or collection-related contamination, including the use of tubes that are not certified to be trace element free. If an elevated result is suspected to be due to contamination, confirmation with a second specimen collected in a certified trace element-free tube is recommended.

Methodology: Inductively Coupled Plasma-Mass Spectrometry (ICP-MS).

Reference Interval:

Test Number	Components	Reference Interval
Arsenic, Whole Blood	Less than or equal to 12.0 $\mu\text{g}/\text{L}$	

Inserted Cells

Inserted Cells