

TEST CHANGE

**Pancreatobiliary FISH
2002528, PF REQUEST**

Specimen Requirements:	
Patient Preparation:	
Collect:	Bile or pancreatic duct brushings, biliary stent, or fine needle aspirates of the pancreas in UroCyteUroVysion-FISH Collection Kit (ARUP Supply #41440) available online through eSupply using ARUP Connect(TM) or contact Client Services at (800) 522-2787. For specific instructions refer to Specimen Collection & Handling.
Specimen Preparation:	Place specimen in Cytolyt or PreservCyt fixative vial. If the specimen is a brushing, submit the brush in the fixative. Specimen should be placed in fixative vial immediately after collection.
Transport Temperature:	Ambient or Refrigerated.
Unacceptable Conditions:	Frozen specimens. Unfixed specimens not collected in Cytolyt or PreservCyt fixative. Specimens submitted in expired collection vials.
Remarks:	Specimen source is required.
Stability:	Ambient: 1 week from collection; Refrigerated: 1 week from collection; Frozen: Unacceptable
Methodology:	Fluorescence in situ Hybridization (FISH)/Computer Assisted Analysis/Microscopy
Performed:	Mon-Fri
Reported:	4-12 days

Note: ~~The UroVysion Kit (ARUP Supply #41440) is designed to detect aneuploidy for chromosomes 3, 7 and 17 via fluorescence in situ hybridization (FISH). Results of this test must be interpreted in conjunction with clinical evidence and other laboratory testing and should not be used alone, as a diagnosis of pancreatobiliary carcinoma.~~
[A positive fluorescence in situ hybridization \(FISH\) result does not identify location or type of malignancy. FISH abnormalities may be associated with high-grade dysplasia or carcinoma in situ. Cytology and biopsy may help clarify such](#)

situations.

CPT Codes: 88366

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:
Refer to Report.

Reference Interval:

<p><u>Negative: No evidence of numeric chromosomal aberrations associated with pancreatobiliary carcinoma identified.</u></p> <p><u>Positive: Numeric chromosomal aberrations associated with pancreatobiliary carcinoma identified.</u></p>	Positive	(≥)5 cells with gains of two or more chromosomes are present
	Equivocal Tetrasomy	Cells present with tetrasomic or near tetrasomic signal patterns
	Equivocal Trisomy	(≥)10 cells with gains of a single chromosome; (Trisomy 7 or Trisomy 3)
	Negative	< 5 abnormal cells are present

Deleted Cells