

HOTLINE: Effective August 19, 2019

2002918	Carbohydrate Deficient Transferrin for Congenital Disorders of Glycosylation (CDG)	CARBOH-CDG
---------	--	------------

Performed: Varies
Reported: 5-13 days

Specimen Required: Collect: Plain red or serum separator tube (SST).
Specimen Preparation: Transfer 0.1 mL serum to an ARUP Standard Transport Tube. (Min: 0.05 mL)
Storage/Transport Temperature: Frozen. Also acceptable: Room temperature or refrigerated.
Remarks: Patient age is required on the test request form. Provide reason (eg, diagnosis) for referral with each specimen.
Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 28 days; Frozen: 45 days