

TEST CHANGE

Alpha-Galactosidase, Serum

2003204, A GALACTO

Specimen Requirements:

Patient Preparation:

Collect: Plain ~~r~~Red or serum separator tube~~Serum Separator Tube~~ (SST).

Specimen Preparation: Transfer 2 mL serum to an ARUP standard transport tube~~Standard Transport Tube~~. (Min: 0.2 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: CRITICAL FROZEN.

Unacceptable Conditions: Thawed specimens.

Remarks: Physician name and phone number are required.

Stability: Ambient: Unacceptable; Refrigerated: 24 hours; Frozen: 2 weeks

Methodology: Quantitative Fluorometry

Performed: Varies

Reported: 8-18 days

Note: Results for this assay are not useful for carrier determination. Carriers usually have levels in the normal range.

CPT Codes: 82657

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report