

TEST CHANGE

Adenovirus, Quantitative PCR	
2007192, ADENO QNT	
Specimen Requirements:	
Patient Preparation:	
Collect:	Lavender (EDTA), pink (K2EDTA), or serum separator tube
Specimen Preparation:	Do not freeze whole blood specimens. Transport 1 mL whole blood, serum or plasma in a sterile container. (Min: 0.5 mL).
Transport Temperature:	Refrigerated
Unacceptable Conditions:	Frozen whole blood. Heparinized specimens.
Remarks:	Specimen source required.
Stability:	Ambient: 24 hours; Refrigerated: 5 days; Frozen: 1 year
Methodology:	Quantitative Real-Time Polymerase Chain Reaction
Performed:	Sun-Sat
Reported:	1-4 days
Note:	The limit of quantification for this DNA assay is 2.3.0 log copies/mL ($2001,000$ copies/mL). If the assay DID NOT DETECT the virus, the test result will be reported as "< 2.3.0 log copies/mL (< $200-1,000$ copies/mL)." If the assay DETECTED the presence of the virus but was not able to accurately quantify the number of copies, the test result will be reported as "Not Quantified."
CPT Codes:	87799
New York DOH Approval Status:	This test is New York DOH approved.
Interpretive Data:	
The quantitative range of this assay is <u>2.</u> 3 <u>-6.3</u> .0- 7.0 log copies/mL (<u>200-2</u> 1,000-10,000,000 copies/mL).	
A negative result (less than 2.3-0 log copies/mL or less than 2001-000 copies/mL) does not rule	

A negative result (less than 2.3.0 log copies/mL or less than 2001,000 copies/mL) does not rule out the presence of PCR inhibitors in the patient specimen or aAdenovirus DNA concentrations below the level of detection of the assay. Inhibition may also lead to underestimation of viral quantitation.



No international standard is currently available for calibration of this assay. Caution should be taken when interpreting results generated by different assay methodologies.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

Reference Interval:

Not detected