

~

TEST CHANGE

.

Leflunomide Metabolite, Serur	n or Plasma
2007460, LEFLUMETSP	
Specimen Requirements:	
Patient Preparation:	Timing of specimen collection: Predose (trough). Obtain specimen 12 - 24 hours after last dose.
Collect:	Plain red. Also acceptable: Lavender (EDTA), pink (K2EDTA), green (sodium heparin), or gray (sodium fluoride).
Specimen Preparation:	Separate from cells within 2 hours of draw. Transfer 1 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.2 mL)
Transport Temperature:	Refrigerated. Also acceptable: Room temperature or frozen.
Unacceptable Conditions:	Whole blood. Potassium oxalate or separator tubes.
Remarks:	
Stability:	Ambient: 7 days; Refrigerated: 17 days; Frozen: 90 days
Methodology:	High Performance Liquid Chromatography-Tandem Mass Spectrometry
Performed:	Sun, Wed, Fri
Reported:	1- <u>7</u> 6 days
Note:	
CPT Codes:	80193
New York DOH Approval Status:	This test is New York DOH approved.

Interpretive Data:

Therapeutic and toxic ranges are not well established. Concentrations greater than 40.000 μ g/mL tend to correlate with improved patient outcome. A proposed therapeutic range is 50.000 - 100.000 μ g/mL. Adverse reactions to Leflunomide, such as diarrhea, hypertension, and liver toxicity, do not correlate well with serum drug concentrations. Leflunomide has a potential risk for teratogenesis. For women being treated with Leflunomide who desire to become pregnant, enhanced drug elimination should be performed until plasma teriflunomide concentrations are lower than 0.020 μ g/mL on two separate tests taken at least 14 days apart.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.



Reference Interval:

Therapeutic	Greater than
Range	40.000 ug/mL
Toxic Level	Not well established.