

TEST CHANGE

Iodine, Serum

2007463, IODINESER

Specimen Requirements:

Patient Preparation: Diet, medication, and nutritional supplements may introduce interfering substances. Patients should be encouraged to discontinue nutritional supplements, vitamins, minerals, nonessential over-the-counter medications for 48 hours (upon the advice of their physician). In addition, the administration of iodine-based contrast media and drugs containing iodine may yield elevated results. During venipuncture, do not use disinfectants (such as Betadine) that contain iodine.

Collect: Royal blue (no additive).

Specimen Preparation: Separate ~~serum~~ from cells ASAP or within 2 hours of collection. ~~Transport~~ 2 mL serum ~~to~~ in an ARUP Trace Element-Free Transport Tube (ARUP supply #43116) available online through eSupply using ARUP Connect(TM) or contact ARUP Client Services at (800-) 522-2787. (Min: 0.5 mL). Do not use utensils (i.e., syringes, needles, or pipettes) in the collection or transfer of the sample. pour directly into transport tube.

Transport Temperature: Refrigerated.

Unacceptable Conditions: Specimens not received in ~~trace element-free~~ ~~Trace Element Free~~ transport tubes. Separator tubes and specimens that are not separated from the clot within 2 hours. Serum collected within 48 hours after administration of a gadolinium (Gd) or iodine (I) containing contrast media (may occur with MRI studies). Plasma.

Remarks:

Stability: Ambient: 30 days; Refrigerated: 30 days; Frozen: 30 days

Methodology: Quantitative Inductively Coupled Plasma-Mass Spectrometry

Performed: Tue, Thu, Sat

Reported: 1-5 days

Note:

CPT Codes: 83018

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Values greater than 250 µg/L may indicate iodine overload.

Elevated results may be due to skin or collection-related contamination, including the use of a noncertified metal-free collection/transport tube. If contamination concerns exist due to elevated levels of serum iodine, confirmation with a second specimen collected in a certified metal-free tube is recommended.

Reference Interval:

40-92 µg/L
