

TEST CHANGE

Adenovirus by Qualitative PCR

2007473, ADENOPCR

Specimen Requirements:

Patient Preparation:

Collect: Lavender (EDTA), pink (K2EDTA), serum separator tube, or urine. Also acceptable: Bronchoalveolar lavage (BAL), CSF, nasopharyngeal swab, sputum, or tissue.

Specimen Preparation: Do not freeze whole blood specimens. Transfer 1 mL whole blood, serum, plasma, BAL, CSF, ~~sputum~~, or urine to a sterile container. (Min: 0.5 mL) Swabs: Transfer to viral transport media ~~2~~ (ARUP supply #12884). Available online through eSupply using ARUP Connect(TM) or contact ARUP Client Services at 800-522-2787. Tissue: Transfer to a sterile container and freeze immediately.

Transport Temperature: Whole blood: Refrigerated. All others: Frozen.

Unacceptable Conditions: Heparinized specimens, tissues in optimal cutting temperature compound.

Remarks: Specimen source required.

Stability: Tissue: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 3 months Urine: Ambient: 3 days; Refrigerated: 14 days; Frozen: 14 days. All others: Ambient: 24 hours; Refrigerated: 5 days; Frozen: 1 year

Methodology: Qualitative Real-Time Polymerase Chain Reaction

Performed: Sun-Sat

Reported: 1-4 days

Note:

CPT Codes: 87798

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test	Components	Reference Interval
------	------------	--------------------

Number		
--------	--	--