

TEST CHANGE

18-Hydroxycorticosterone by Mass Spectrometry

2007582, 18 HYDRO

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube (SST) or plasma preparation tube (PPT).
Also acceptable: Lavender (EDTA) or green (sodium heparin).

Specimen Preparation: Separate from cells within 1 hour of collection. Transfer 3 mL serum or plasma to an ARUP standard transport tube and freeze immediately. (Min: 1 mL)
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: CRITICAL FROZEN.

Unacceptable Conditions:

Remarks:

Stability: Ambient: Unacceptable24 hours; Refrigerated: Unacceptable24 hours; Frozen: 3 months

Methodology: Quantitative Mass Spectrometry

Note:

CPT Codes: 82542

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Refer toBy report