

**TEST CHANGE**

**Encephalitis Panel with Reflex to Herpes Simplex Virus Types 1 and 2 Glycoprotein G-Specific Antibodies, IgG, Serum**

2008915, ENCEPH

**Specimen Requirements:**

**Patient Preparation:**

**Collect:** Serum separator tube.

**Specimen Preparation:** Transfer ~~4.0mL~~~~3-mL~~ serum to an ARUP [standard transport tube](#). ~~Standard Transport Tube~~. (Min: ~~2.0mL~~~~1.5 mL~~)

**Transport Temperature:** Refrigerated.

**Unacceptable Conditions:** Refer to individual components. CSF (refer to Encephalitis Panel with Reflex to Herpes Simplex Virus Types 1 and 2 Glycoprotein G-Specific Antibodies, IgG, CSF, test code 2008916).

**Remarks:**

**Stability:** After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year

**Methodology:** Semi-Quantitative Enzyme-Linked Immunosorbent Assay ([ELISA](#))/Semi-Quantitative Chemiluminescent Immunoassay ([CLIA](#))

**Performed:** Sun-Sat

**Reported:** 2-6 days

**Note:** If HSV 1 and/or 2 IgG is 1.10 IV or greater, then HSV 1 G-Specific IgG and HSV 2 G-Specific IgG will be added. Additional charges apply.

**CPT Codes:** 86765 x2; 86735 x2; 86787 x2; 86789; 86788; 86694 x2; if reflexed, add 86695; 86696

**New York DOH Approval Status:** This test is New York DOH approved.

**Interpretive Data:**

Component	Interpretation
Measles (Rubeola) Antibody, IgG	13.4 AU/mL or less: Negative - No significant level of detectable measles (rubeola) IgG antibody. 13.5-16.4 AU/mL: Equivocal - Repeat testing in 10-14 days may be helpful. 16.5 AU/mL or greater: Positive - IgG antibody to measles (rubeola) detected, which may indicate a current or past exposure/immunization to measles (rubeola).
Measles (Rubeola) Antibody, IgM	0.79 AU or less: Negative - No significant level of IgM antibodies to measles (rubeola) virus detected. 0.80-1.20 AU: Equivocal - Repeat testing in 10-14 days may be helpful. 1.21 AU or greater: Positive - IgM antibodies to measles (rubeola) virus detected. Suggestive of current or recent infection or immunization. However, low levels of IgM antibodies may occasionally persist for more than 12 months post-infection or immunization.
Mumps Virus Antibody, IgG	8.9 AU/mL or less: Negative - No significant level of detectable IgG mumps virus antibody. 9.0-10.9 AU/mL: Equivocal - Repeat testing in 10-14 days may be helpful. 11.0 AU/mL or greater: Positive - IgG antibody to mumps virus detected, which may indicate a current or past exposure/immunization to mumps virus.
Mumps Virus Antibody, IgM	0.79 IU or less: Negative - No significant level of detectable IgM

	<p>antibody to Mumps virus. 0.80-1.20 IV: Equivocal - Borderline levels of IgM antibody to Mumps virus. Repeat testing in 10-14 days may be helpful. 1.21 IV or greater: Positive - Presence of IgM antibody to Mumps virus detected, which may indicate a current or recent infection. However, low levels of IgM antibody may occasionally persist for more than 12 months post-infection or immunization.</p>
Varicella-Zoster Virus Antibody, IgG	<p>134.9 IV or less: Negative - No significant level of detectable varicella-zoster IgG antibody. 135.0-164.9 IV: Equivocal - Repeat testing in 10-14 days may be helpful. 165.0 IV or greater: Positive - IgG antibody to varicella-zoster detected, which may indicate a current or past varicella-zoster infection.</p>
Varicella-Zoster Virus Antibody, IgM	<p>0.90 ISR or less: Negative - No significant level of detectable varicella-zoster virus IgM antibody. 0.91-1.09 ISR: Equivocal - Repeat testing in 10-14 days may be helpful. 1.10 ISR or greater: Positive - Significant level of detectable varicella-zoster virus IgM antibody. Indicative of current or recent infection. However, low levels of IgM antibodies may occasionally persist for more than 12 months post-infection or immunization.</p>
Herpes Simplex Virus Type 1	<p>0.89 IV or less: Not Detected. 0.90-1.09 IV: Indeterminate - Repeat</p>

and/or 2 Antibodies, IgM by ELISA	testing in 10-14 days may be helpful. 1.10 IV or greater: Detected - IgM antibody to HSV detected, which may indicate a current or recent infection. However, low levels of IgM antibodies may occasionally persist for more than 12 months post-infection.
Herpes Simplex Virus Type 1 and/or 2 Antibodies, IgG	0.89 IV or less: Not Detected. 0.90-1.09 IV: Indeterminate - Repeat testing in 10-14 days may be helpful. 1.10 IV or greater: Detected.
West Nile Virus Antibody, IgG by ELISA, Serum	1.29 IV or less: Negative - No significant level of West Nile virus IgG antibody detected. 1.30-1.49 IV: Equivocal - Questionable presence of West Nile virus IgG antibody detected. Repeat testing in 10-14 days may be helpful. 1.50 IV or greater: Positive - Presence of IgG antibody to West Nile virus detected, suggestive of current or past infection.
West Nile Virus Antibody, IgM by ELISA, Serum	0.89 IV or less: Negative - No significant level of West Nile virus IgM antibody detected. 0.90-1.10 IV: Equivocal - Questionable presence of West Nile virus IgM antibody detected. Repeat testing in 10-14 days may be helpful. 1.11 IV or greater: Positive - Presence of IgM antibody to West Nile virus detected, suggestive of current or recent infection.

**Reference Interval:**

Test Number	Components	Reference Interval
	HSV 1 and/or 2 Abs, IgM by ELISA	0.89 IV or less
	West Nile Virus Ab, IgG, Ser	1.29 IV or less
	West Nile Virus Ab, IgM, Ser	0.89 IV or less
	Mumps Virus Antibody, IgM	0.79 IV or less
	Measles, Rubella, Antibody IgM	0.79 AU or less
	Varicella-Zoster Virus Antibody, IgM	0.90 ISR or less