

Quarterly HOTLINE: Effective February 19, 2019

2009373 Zinc Quantitative, Whole Blood ZINC WB

Performed: Sun-Sat **Reported:** 1-3 days

Specimen Required: Patient Prep: Diet, medication, and nutritional supplements may introduce interfering substances. Patient should be encouraged to

discontinue nutritional supplements, vitamins, minerals, and non-essential over-the-counter medications (upon the advice of their

physician).

Collect: Royal Blue (K₂EDTA or Na₂EDTA).

Specimen Preparation: Transport 2 mL whole blood in the original collection tube. (Min: 0.5 mL)

Storage/Transport Temperature: Room temperature. Also acceptable: Refrigerated.

Unacceptable Conditions: Specimens collected in tubes other than Royal Blue (EDTA). Specimens transported in containers other

than a Royal Blue (EDTA) tube or Trace Element-Free Transport Tube. Heparin anticoagulant. Clotted specimens. Stability (collection to initiation of testing): Ambient: Indefinitely; Refrigerated: Indefinitely; Frozen: Unacceptable

Reference Interval: Effective February 19, 2019

 $440.0-860.0 \, \mu g/dL$

Interpretive Data: Elevated results may be due to skin or collection-related contamination, including the use of a noncertified metal-free collection/transport tube. If contamination concerns exist due to elevated levels of blood zinc, confirmation with a second specimen collected in a certified metal-free tube is recommended.

Zinc concentration in blood has not been shown to change significantly in deficiency or with supplementation.

See Compliance Statement B: www.aruplab.com/CS

HOTLINE NOTE: There is a numeric map change associated with this test.

Change the numeric map for component 2009374, Zinc Quantitative, Whole Blood from XXXXXX.X to XXXX.X.