

TEST CHANGE

Histoplasma Galactomannan Antigen Quantitative by EIA, Urine

2009418, HISTOGM U

Specimen Requirements:

Patient Preparation:

Collect: Random urine.

Specimen Preparation: Transfer 2 mL urine to an ARUP [Standard Transport Tube](#) ~~standard transport tube~~.

Transport Temperature: Refrigerated.

Unacceptable Conditions: Specimens other than urine. Urine in boric acid. Serum; refer to test Histoplasma Antigen by EIA, Serum (ARUP test code 0092522).

Remarks:

Stability: Ambient: 24 hours; Refrigerated: 2 weeks; Frozen: 2 weeks

Methodology: Quantitative Enzyme Immunoassay [\(EIA\)](#)

Performed: Sun-Sat

Reported: 1-2 days

Note:

CPT Codes: 87385

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Less than 0.4 ng/mL ~~??~~ = Not Detected

0.4-0.7 ng/mL ~~??~~ = Detected (below the limit of quantification)

0.8-24.0 ng/mL ~~??~~ = Detected

Greater than 24.0 ng/mL ~~??~~ = Detected (above the limit of quantification)

The quantitative range of this assay is 0.8-24.0 ng/mL. Antigen concentrations between 0.4-0.7 or ~~??~~ > 24.0 ng/mL fall outside the linear range of the assay and cannot be accurately quantified.

This EIA test should be used in conjunction with other diagnostic procedures, including microbiological culture, histological examination of biopsy samples, and/or radiographic evidence, to aid in the diagnosis of histoplasmosis.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

Reference Interval:

Not Detected

HOTLINE NOTE: There is a result type change associated with this test. The result type has changed from non-numeric to numeric or vice versa. Refer to the Hotline Test Mix for interface build information.

HOTLINE NOTE: There is a numeric map change associated with this test. Refer to the Hotline Test Mix for interface build information.