

TEST CHANGE

Alcohol, Urine, Quantitative

2010136, CDCO ETOH

Specimen Requirements:

Patient Preparation:

Collect: Random urine.

Specimen Preparation: Transfer 4 mL urine without additives or preservatives to an ARUP Standard Transport Tube. (Min: 1 mL)

Transport Temperature: Room temperature.

Unacceptable Conditions:

Remarks:

Stability: Ambient: 1 week; Refrigerated: 1 month; Frozen: 3 years (Avoid repeated freeze/thaw cycles)

Methodology: Quantitative Gas Chromatography

Performed: Sun-Sat

Reported: 1-4 days

Note:

CPT Codes: 80320 (Alt code: G0480)

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Methodology: Quantitative Gas Chromatography with Flame Ionization Detection

Positive cutoff: 5 mg/dL

For medical purposes only; not valid for forensic use.

The absence of expected drug(s) and/or drug metabolite(s) may indicate inappropriate timing of specimen collection relative to drug administration, poor drug absorption, diluted/adulterated urine, or limitations of testing. The concentration value must be greater than or equal to the cutoff to be reported as positive. Interpretive questions should be directed to the laboratory.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

Reference Interval:

Effective August 17, 2015

Drugs Covered	Cutoff Concentrations	
Ethanol	5 mg/dL	
Ethyl Alcohol	5 mg/dL	

HOTLINE NOTE: There is a numeric map change associated with this test. Refer to the Hotline Test Mix for interface build information.