

TEST CHANGE

Thrombomodulin by Immunohistochemistry

2010170, THROM IHC

Specimen Requirements:

Patient Preparation:

Collect: Tissue or cells.

Specimen Preparation: Formalin fix (10 percent neutral buffered formalin) and paraffin embed specimen (cells must be prepared into a cellblock). Protect paraffin block and/or slides from excessive heat. Transport tissue block or 3 unstained (4-micron thick sections), positively charged slides in a tissue transport kit (ARUP supply #47808) available online through eSupply using ARUP Connect(TM) or contact ARUP Client Services at (800-) 522-2787. (Min: 3 slides) If sending precut slides, do not oven bake. Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Room temperature. Also acceptable: Refrigerated. Ship in cooled container during summer months.

Unacceptable Conditions: Specimens submitted with ~~nonrepresentative~~
~~representative~~ tissue type. Depleted specimens.

Remarks:

Stability: Ambient: Indefinitely; Refrigerated: Indefinitely; Frozen: Unacceptable

Methodology: Immunohistochemistry (IHC)

Performed: Varies

Reported: ~~3-4~~7 days

Note: All stains will be handled as "Stain and Return" without interpretation.

CPT Codes: 88342

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval
-------------	------------	--------------------