

HOTLINE: Effective February 1, 2021

## 2011148 Herpes Simplex Virus (HSV) by PCR with Reflex to HSV (HSV-1/HSV-2) HSVPCR RFX Subtype by PCR

Specimen Required: Collect: Lavender (EDTA), pink (K2EDTA), serum separator tube. OR CSF, bronchoalveolar lavage (BAL), amniotic fluid, vesicle

fluid, ocular fluid, tissue. **OR** endocervical specimen in ThinPrep Pap Test media.

Specimen Preparation: Separate plasma or serum from cells. Transfer 1 mL plasma, serum, CSF, BAL, amniotic fluid, ocular fluid or ThinPrep specimen to a sterile container. (Min: 0.5 mL)

**Tissue:** Transfer to a sterile container and freeze immediately.

Vesicle Fluid: Transfer to viral transport media (ARUP supply #12884). Available online through eSupply using ARUP Connect™ or contact ARUP Client Services at (800) 522-2787.

<u>Storage/Transport Temperature:</u> Frozen. <u>Remarks:</u> **Specimen source required.** 

<u>Unacceptable Conditions:</u> Heparinized specimens, tissues in optimal cutting temperature compound.

Stability (collection to initiation of testing): Tissue: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 3 months

All others: Ambient: 8 hours; Refrigerated: 72 hours; Frozen: 3 months