

Quarterly HOTLINE: Effective February 19, 2019

2011521	Meprobamate, Serum or Plasma, Quantitative	MEPRO SP
----------------	---	-----------------

Performed: Sun, Tue, Thu
Reported: 1-4 days

Specimen Required: Collect: Plain Red, Lavender (K₂EDTA), Lavender (K₃EDTA), or Pink (K₂EDTA).
Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.5 mL)
Storage/Transport Temperature: Refrigerated.
Unacceptable Conditions: Whole blood. Gel Separator Tubes, Light Blue (Sodium Citrate), or Yellow (SPS or ACD Solution).
Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: 3 months

Interpretive Data: The therapeutic range is based on serum pre-dose (trough) draw at steady-state concentration. Adverse effects may include drowsiness, dizziness and headache.

See Compliance Statement B: www.aruplab.com/CS

HOTLINE NOTE: Remove information found in the Specimen Required Remarks field. There is also a component change associated with this test.

Remove component 2011522, Meprobamate Dose
 Remove component 2011524, Meprobamate Route
 Remove component 2011523, Meprobamate Dose Frequency
 Remove component 2011525, Meprobamate Type of Draw