

## TEST CHANGE

### THC (Cannabinoids), Urine Screen with Reflex to Quantitation

2012270, THC RFX U

#### Specimen Requirements:

##### Patient Preparation:

Collect: Random urine.

Specimen Preparation: Transfer 4 mL urine with no additives or preservatives to an ARUP Standard Transport Tube. (Min: 2 mL)

Transport Temperature: Room temperature.

Unacceptable Conditions: Breast milk.

##### Remarks:

Stability: Ambient: 1 week; Refrigerated: 1 month; Frozen: 1 month

Methodology: Qualitative Enzyme Immunoassay/Quantitative Liquid Chromatography-Tandem Mass Spectrometry

Performed: Sun-Sat

Reported: 1-4 days

Note: If the specimen screens positive, then Confirmation/Quantitation by LC-MS/MS (ARUP test code 0090369) will be added to confirm result. Additional charges apply. Compare to Pain Management, Marijuana Metabolite, Quantitative, with medMATCH, Urine; Pain Management, Marijuana Metabolite, with Confirmation with medMATCH, Urine.

CPT Codes: 80307; if reflexed, add 80349 (Reflexed Alt Code: G0480)

New York DOH Approval Status: This test is New York DOH approved.

#### Interpretive Data:

The absence of expected drug(s) and/or drug metabolite(s) may indicate non-compliance, inappropriate timing of specimen collection relative to drug administration, poor drug absorption, diluted/adulterated urine, or limitations of testing. The concentration at which the screening test can detect a drug or metabolite varies. Specimens for which drugs or drug classes are detected by the screen are reflexed to a second, more specific technology (GC/MS and/or LC-MS/MS). The concentration value must be greater than or equal to the cutoff to be reported as positive. Interpretive questions should be directed to the laboratory.

For medical purposes only; not valid for forensic use.

Reference Interval:

Screen cutoff concentration: ~~50~~<sup>20</sup> ng/mL

---