

HOTLINE: Effective May 17, 2021

2012618

Risk of Ovarian Malignancy Algorithm

ROMA

Specimen Required: Collect: Serum Separator Tube (SST). Also acceptable: Green (Sodium or Lithium Heparin), Lavender (EDTA), or Pink (K₂ EDTA).
Specimen Preparation: Allow specimen to clot completely at room temperature. Transfer 1.5 mL serum to an ARUP Standard Transport Tube. (Min: 1 mL)
Storage/Transport Temperature: Frozen.
Unacceptable Conditions: Hemolyzed specimens.
Stability (collection to initiation of testing): Ambient: 5 hours; Refrigerated: 48 hours; Frozen: 3 months