

Quarterly HOTLINE: Effective November 12, 2018

**2012674 Human Immunodeficiency Virus (HIV) Combo Antigen/Antibody (HIV-1/O/2) by HIV PANEL
CIA, Reflexive Panel**

Methodology: Qualitative Chemiluminescent Immunoassay/Qualitative Immunoassay/Quantitative **Transcription-Mediated Amplification**
Performed: Sun-Sat
Reported: 1-2 days

Specimen Required: Collect: Lavender (EDTA) or Pink (K₂EDTA).

Specimen Preparation: Separate from **cells within 24** hours of collection. Transfer 3 mL plasma into an ARUP Standard Transport Tube. (Min: 2 mL) Remove particulate material. This test requires a dedicated transport tube submitted only for HIV testing.

Storage/Transport Temperature: Frozen.

Unacceptable Conditions: Serum. Heparinized or citrated plasma specimens. Plasma preparation tube. Specimens containing particulate material. Severely hemolyzed or heat-inactivated specimens.

Stability (collection to initiation of testing): After separation from cells: Ambient: Unacceptable; Refrigerated: **72 hours**; Frozen: **3 months** (avoid repeated freeze/thaw cycles)

Reference Interval:

Effective November 12, 2018

Test Number	Components	Reference Interval		
	HIV 1,2 Combo Antigen/Antibody	Negative		
2012669	Human Immunodeficiency Virus Types 1 and 2 (HIV-1/2) Antibody Differentiation, Supplemental, with Reflex to HIV-1 Quantitative NAAT, Plasma	Test Number	Components	Reference Interval
			HIV-1 Antibody	Negative
			HIV-2 Antibody	Negative
		3000867	Human Immunodeficiency Virus 1 (HIV-1) by Quantitative NAAT, Plasma	Not detected

Note: The fourth-generation screen test is for the simultaneous qualitative detection of Human Immunodeficiency Virus Type 1 (HIV-1) p24 antigen and antibodies to HIV Type 1 (HIV-1 groups M and O) and HIV Type 2 (HIV-2). Results of the screen cannot be used to distinguish between the presence of HIV-1 p24 antigen, HIV-1 antibody, or HIV-2 antibody.

If the HIV-1,2 Combo Antigen/Antibody screen is repeatedly reactive, then the HIV-1/2 Ab Differentiation Immunoassay will be performed. Additional charges apply. The HIV-1/2 Ab Differentiation Immunoassay confirms and discriminates between HIV-1 and HIV-2 antibodies. Results for each type are reported.

If the HIV-1/2 Ab Differentiation Immunoassay is Negative or Indeterminate, then the Human Immunodeficiency Virus 1 (**HIV-1**) by Quantitative **NAAT, Plasma** will be added. Additional charges apply.

This multi-test algorithm is recommended by the Centers for Disease Control and Prevention (CDC) and was adopted by the Clinical Laboratory Standards Institute (CLSI) for the diagnosis of HIV (refer to <http://www.arupconsult.com/Topics/HIV.html>).

Refer to the following tests for additional information regarding Performed or Reported times, Interpretive Data and Notes for the reflex tests of this panel: Human Immunodeficiency Virus Types 1 and 2 (HIV-1/2) Antibody Differentiation, Supplemental with Reflex to HIV-1 Quantitative **NAAT, Plasma** (ARUP test code 2012669); Human Immunodeficiency Virus 1 (**HIV-1**) by Quantitative **NAAT, Plasma** (ARUP test code 3000867).

HOTLINE NOTE: There is a reflexive pattern change associated with this test.

Add reflex to 3000867, HIV-1 by Quantitative NAAT, Plasma to reflexive orderable 2012669

Remove reflex to 0055598, HIV-1 by Quantitative PCR from reflexive orderable 2012669