

HOTLINE: Effective May 17, 2021

2013662 Cystic Fibrosis (CFTR) 165 Pathogenic Variants, Fetal

CF VAR FE

Specimen Required: Collect: Fetal Specimen: Two T-25 flasks of cultured amniocytes at 80 percent confluency. *If the client is unable to culture

amniocytes, this can be arranged by contacting ARUP Client Services at (800) 522-2787.

Maternal Specimen: Lavender (EDTA), Pink (K2EDTA), or Yellow (ACD Solution).

Specimen Preparation: Cultured Amniocytes: Fill flasks with culture media. Backup cultures must be retained at the client's

institution until testing is complete.

Maternal Specimen: Transport 3 mL whole blood. (Min. 1 mL)

Storage/Transport Temperature: CulturedAmniocytes: CRITICAL ROOM TEMPERATURE. Must be received within 48 hours

of shipment due to lability of cells.

MaternalSpecimen: Refrigerated.

Remarks: Maternal sample is recommended for proper test interpretation; order Maternal Cell Contamination, Maternal

Specimen. Patient History Form is available on the ARUP Web site or by contacting ARUP Client Services.

<u>Unacceptable Conditions:</u> Maternal Specimen: Plasma or serum. Specimens collected in sodium heparin or lithium heparin tubes.

Frozen specimens in glass collection tubes.

Stability (collection to initiation of testing): Fetal Specimen: Ambient: 48 hours; Refrigerated: Unacceptable; Frozen: Unacceptable

Maternal Specimen: Ambient: 72 hours; Refrigerated: 1 week; Frozen: 1 month