

HOTLINE: Effective June 3, 2019

3000202 **5-Hydroxyindoleacetic acid (5-HIAA), Plasma****5 HIAA PLA****Performed:** Varies
Reported: 14-17 days**Specimen Required:** Patient Prep: Fast overnight prior to collection.Collect: Z plasma preservative tube (ARUP supply #40874) available online through eSupply using ARUP Connect™ or contact ARUP Client Services at (800) 522-2787.Specimen Preparation: Separate from cells within 10 minutes. Transfer 4 mL plasma to an ARUP Standard Transport Tube and freeze immediately. (Min: 1 mL)Storage/Transport Temperature: CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.Unacceptable Conditions: Thawed specimens.Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 1 week; Frozen: 6 months