

Quarterly HOTLINE: Effective August 20, 2018

New Test **3000576** **Hepatitis C Virus (HCV) by Quantitative NAAT with Reflex to HCV Genotype by Sequencing** **HCVQT GR**



Additional Technical Information

Methodology: Quantitative Transcription Mediated Amplification/Sequencing
Performed: Sun-Sat
Reported: 1-8 days

Specimen Required: Collect: Lavender (EDTA), Pink (K₂EDTA), Plasma Preparation Tube (PPT), or Serum Separator Tube (SST).
Specimen Preparation: Separate from cells within 6 hours of collection. Transfer 3 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 1.7 mL)
Storage/Transport Temperature: Frozen.
Unacceptable Conditions: Heparinized specimens.
Stability (collection to initiation of testing): After separation from cells: Ambient: Unacceptable; Refrigerated: 72 hours; Frozen: 6 weeks

Reference Interval:

Available Separately	Components	Reference Interval
3000572	Hepatitis C Virus by Quantitative NAAT	Not Detected
0055593	Hepatitis C Virus Genotype by Sequencing	By report

Interpretive Data: Refer to report.

Note: If Hepatitis C Virus by Quantitative NAAT result is greater than or equal to 4,000 IU/mL, then Hepatitis C Virus Genotype by Sequencing will be added. Additional charges apply.

CPT Code(s): 87522; if reflexed, add 87902

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.