

HOTLINE: Effective August 17, 2020

New Test	3001238	Human Immunodeficiency Virus Type 1 (HIV-1) GenoSure PRIme	HIV GSPRI
-----------------	----------------	---	------------------

Methodology: Polymerase Chain Reaction/Sequencing
Performed: Varies
Reported: 7-13 days

Specimen Required: Collect: Lavender (EDTA) or Plasma Preparation Tube (PPT).
Specimen Preparation: Separate from cells within 6 hours of collection. Transfer 5 mL plasma to ARUP Standard Transport Tubes and freeze immediately. (Min: 3 mL)
Storage/Transport Temperature: **CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.**
Unacceptable Conditions: Thawed specimens.
Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 2 years

Reference Interval: By Report

Note: Procedure should be used for patients with documented HIV-1 infection and viral loads greater than 500 copies/mL.

CPT Code(s): 87900; 87901; 87906

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.